

# Aquagent®

The Scharlau comprehensive pyridine-free solutions range for a reliable Karl Fischer Titration





Karl Fischer titration is a well known and globally accepted method for water determination since the beginning of the 20th century. KF titration uses volumetric or coulometric titration to determine the water content in a wide variety of samples including chemicals, pharmaceuticals, food and oils. It is used both in industrial processes as well as in quality control laboratories.

The first KF reagents that were developed contained pyridine, which was assumed to be essential for the reaction. Further experiments demonstrated that pyridine could be replaced by other basic compounds, which were able to play the same role providing less toxicity.

Most of the pyridine-free reagents, including our Aquagent®, contain imidazole instead of pyridine. Imidazole is a non-toxic base, has a good buffering capacity and allows fast and stable titration end-points.



# Aquagent®: An appropriate product range

AQUAGENT® is the comprehensive product range of Scharlau pyridine-free volumetric and coulometric Karl Fischer reagents for water analysis. We offer a wide range of safe, reliable and easy-to-use pyridine-free solutions and standards to meet any requirement of current laboratories for Karl Fischer titration.

#### Aquagent® comprises:

- Volumetric titration with One and Two-component reagents
- ♠ Coulometric titration with cells with and without diaphragm
- Water Standards

#### Suitable for the following applications:

- Ketones and Aldehydes
- Carbohydrates, Inorganic Salts and Proteins
- Oils and Fats
- Crude and related products
- Strong acids

# Why to use Aquagent®

Water content can affect product quality, texture, shelf life, chemical stability and reactivity. So Aquagent® provides an accurate water content determination in volumetric and coulometric titration with unique performance.

- Less toxicity, more safety
- End point stability
- Accuracy and reproducibility
- Faster titration
- No bad and noxious smells
- Wide range of capacities
- Decreased environmental impact
- Wide applicability
- Outstanding quality: quality control under rigorous standards
- Selected raw materials
- ♦ Globally available: international sales network
- 30 years experience

Aquagent® providing you reliable results in volumetric and coulometric Karl Fischer titration





Two methods are available for the determination of water content based on the Karl Fischer reaction: volumetric and coulometric. The choice of the method primarily depends of the amount of water expected in the sample. It is known that the choice of the right product is a key factor in obtaining reliable and reproducible results.

# Aquagent® volumetric solutions

In case of higher water content (0,1- 100%), the volumetric titration is the method to choice. It is the most used for water titration, which is determined by measuring the required volume of Karl Fischer reagent consumed to reach the titration end-point. This end-point is indicated by an excess of iodine and measured potentiometrically.

Scharlau supplies both one-component and two-component reagents for volumetric titration.



### Aquagent® One-component

In one-component Karl Fischer titration, all substances involved in the reaction are mixed in one reagent: the titrant. The one-component reagents are user-friendly and allow more flexibility in the choice of the more suitable solvent according to the type of sample; on the other hand, they should be frequently titrated due to the reactivity of their components.

Scharlau offers a range of one-component solutions suitable for general use as well as for specific applications.

#### **GENERAL USE**

#### Reagents:

#### **Aquagent® Complet 2**

A general purpose reagent for samples with low and medium water content. Each mL titrates approx. 2 mg water. Generally used together with methanol as a solvent.

#### Aquagent® Complet 5

A general purpose reagent for samples with high and medium water content. Each mL titrates approx. 5 mg water. Generally used together with methanol as a solvent.

	DESCRIPTION	CAPACITY	ART. NO.
		500mL	AQ00070500
	Aquagent® Complet 2	1L	AQ00071000
		2,5L	AQ00072500
		500mL	AQ00030500
	Aquagent® Complet 5	1L	AQ00031000
		2.51	AQ00032500

#### Solvents:

#### **Dry Methanol**

The sample should always be dissolved in a dry solvent to be titrated. The most common is dry methanol. If the sample is not soluble in methanol, any other dry solvent should be used (see specific applications on next page).

#### **Aquagent® Methanol Fast**

An improved formulation for a faster KF volumetric titration.





FIGURE 1: Shows the time to achieve the end-point in the KF volumetric one-component titration using different types of methanol solvent. Sample:  $20 \text{ mg H}_2\text{O}$  injection by weight.

DESCRIPTION	CAPACITY	ART. NO.
Methanol, dry (max. 0,005% H₂O),	1L	ME03041000
reagent grade	2,5L	ME03042500
A	1L	AQ00111000
Aquagent® Methanol Fast	2,5L	AQ00112500





# Aquagent® volumetric solutions



Aquagent® One-component

#### SPECIFIC APPLICATIONS

#### Aquagent® Complet 5K

Aldehydes and ketones react with methanol releasing water as by-product of this reaction. Hence, when the sample contains aldehydes or ketones erroneous results might be obtained. To avoid this effect a specific reagent is needed: our Aquagent® Complet 5K. It is used in conjunction with Aquagent® Medium K, a specific solvent that does not contain methanol. The usual titre is 5 mg water/mL.

DESCRIPTION	CAPACITY	ART. NO.
DESCRIPTION  Aquagent® Complet 5K	500mL	AQ00040500
	1L	AQ00041000

#### Aquagent® Medium K

Methanol reacts with both ketones and aldehydes and water is a by-product of these reactions. For this reason, when the sample contains aldehydes or ketones, methanol shall be substituted by another solvent, our Aquagent® Medium K.


#### Aquagent® Buffer Acid (additive)

The Karl-Fischer reaction optimally runs at pH values between 5 and 7. When determining water in strongly acidic compounds, it is recommended to neutralize the working medium with our Aquagent® Buffer Acid.

#### **Dry formamide (additive)**

Formamide improves the solubility of carbohydrates, proteins and inorganic salts. This solvent can be added to methanol in no more than 50% by volume.

DESCRIPTION	CAPACITY	ART. NO.
Aguagent® Buffer, acid	500mL	AQ00090500 AQ00091000 AQ00050500 AQ00051000
Aquagent <sup>a</sup> Buller, acid	1L	AQ00091000
Agus goot® Madium I/	500mL	AQ00050500
Aquagent® Medium K	1L	AQ00051000
Formamide, dry (max. 0,02% H <sub>2</sub> O), reagent grade	1L	FO00281000



# **Aquagent® One-component reagents:**

- ♦ High titration rate for fast analyses
- ♦ Ensure reproducible and consistent high quality results
- Flexibility: the solvent can be suited to the sample matrix
- Unlimited water capacity compared to two components reagent







### Aquagent® Two-component

In two-component systems, the solvent component does not act just as solvent medium, but also contains part of the reagents. This allows longer shelf-life and avoids the need for frequent re-titration.

The use of two component reagents is more expensive but offers advantages, compared to one-component reagents: faster titration, less consumption of titration reagents and better long-term stability of the reagents.

Scharlau offers a range of titrants and solvents which are suitable for general use as well as for specific applications.

#### **GENERAL USE**

#### Aquagent® Titrant 2

A general purpose reagent that contains iodine and methanol. Titre is aprox. 2 mg water/mL. Shall be used in conjunction with Aquagent® Solvent.

#### **Aquagent® Titrant 5**

A general purpose reagent that contains iodine and methanol. Titre is approx. 5 mg water/mL. Shall be used in conjunction with Aquagent® Solvent.

DESCRIPTION	CAPACITY	ART. NO.
Aguagant® Titrant 2	500mL	AQ00060500
Aquagent® Titrant 2	1L	AQ00061000
	500mL	AQ00010500
Aquagent® Titrant 5	1L	AQ00011000
	2,5L	AQ00012500

Scharlau offers several products to be used as the solvent-component in conjunction with Aquagent® Titrant:

#### Aquagent® Solvent

A general reagent that contains SO<sub>2</sub>, imidazole and methanol. Shall be used use in conjunction with Aquagent® Titrant.

# Scharlau

DESCRIPTION	CAPACITY	ART. NO.
A successful Columns	1L	AQ00021000
Aquagent® Solvent	2,5L	AQ00022500

Aquagent® Metanol Fast

#### SPECIFIC APPLICATIONS

#### Aquagent® Solvent CM

Solvent-component for titration of fats and oils, with additives to facilitate solubility of long-chained hydrocarbons.

DESCRIPTION	CAPACITY	ART. NO.
Aguagant® Calvent CM	1L	AQ00081000
Aquagent® Solvent CM	2.5L	AQ00082500

#### Aquagent® Solvent OIL

Solvent-component for titration of fats and oils, free from halogenated hydrocarbons.

DESCRIPTION	CAPACITY	ART. NO.
Aquagent® Solvent OIL	1L	AQ00101000
7 Iquagorii Gorvani Giz		710001010

# **Aquagent® Two-component reagents:**

- **♦** Faster titration in comparison to one-component reagents
- High accuracy for high quality results
- Titre more stable in comparison to one-component reagents





Coulometric Karl Fischer titration is indicated for low water content at ppm level (<0,1%) or for water determination in valuable samples. In coulometric titration, the iodine required for the reaction is formed at the electrode in the titration vessel itself by anodic oxidation. The water content is accurately calculated from the current used over a specific time period. The measuring cell contains an anode and a cathode compartment which can be separated by a membrane or diaphragm. The titrators cells can therefore be with or without diaphragm.

# **Aquagent® Coulometric solutions**

Scharlau offers a suitable AQUAGENT® for both cell types.

# Aquagent® for cells with diaphragm

#### Anolyte:

# Aquagent® Coulometric A Anolyte for coulometric KF titration

Suitable for cells with diaphragm. This general purpose reagent contains the reaction compounds for the anodic compartment from the electrolytic cell. Should be used in conjuntion with Aquagent® Coulometric CG.

# Aquagent® Coulometric Oil Anolyte for coulometric KF titration

Suitable for cells with diaphragm. This reagent for the anodic compartment is specially formulated for crude and related products samples. Should be used in conjunction with Aquagent® Coulometric CG.

#### Catholyte:

# Aquagent® Coulometric CG Catholyte for coulometric KF titration

Suitable for cells with diaphragm. This reagent contains the reaction compounds for the catodic compartment from the electrolytic cell. Should be used in conjunction with Aquagent® Coulometric A or Oil.

DESCRIPTION	CAPACITY	ART. NO.
Aquagent® Coulometric A, anolyte	500mL	AQ00220500
Aquagent® Coulometric Oil, anolyte	100mL	AQ00250100
Aquagent® Coulometric CG,	100mL	AQ00230100
catholyte	10x5mL	AQ00230050



# Aquagent® for cells without diaphragm

# Aquagent® Coulometric AG, for coulometric KF titration

Suitable for cells without diaphragm.

One single solution containg all reaction compounds.

DESCRIPTION	CAPACITY	ART. NO.
Aquagent® Coulometric AG	500mL	AQ00240500
	1L	AQ00241000





## Scharlau Standards for Karl Fischer titration

Standards of a known water content are used to determine the titre of the reagents. They are more and more requested due to an increasing demand for more reliable and comparable results. Our Aquagent® product family includes:

- · Solid standards: sodium tartrate dihydrate, stable, non-hygroscopic, with a water content about 15,66%.
- Liquid standards: Aquagent® Standard 1.0 for coulometric and Aguagent® Standard 10.0 for volumetric Karl Fischer titrations respectively. We pack our 1.0 and 10.0 standards in vials to maintain optimum conditions until they are opened. Each vial provides sufficient standard for one titration. Aquagent® Standard 5.0, suitable for daily titre control as well as for equipment validation. Shelf life is stable and included in the Certificate of Analysis.

DESCRIPTION	CAPACITY	ART. NO.
Aquagant® di Cadium tartrata dibudrata	25g	AQ00300025
Aquagent® di-Sodium tartrate dihydrate	100g	AQ00300100
Aquagent® standard solution 1.0* (1mg/g)	10 x 4mL	AQ00190040
Aquagent® standard solution 10.0* (10mg/g)	10 x 8mL	AQ00200080
Aguagant® standard solution F.O. (Fmg./ml.)	100mL	AQ00210100
Aquagent® standard solution 5.0 (5mg/mL) -	500mL	AQ00210500

<sup>\*</sup>Traceable to NIST



# **Benefits of Aquagent®**

- Highest quality results
- Increased safety
- Time saving
- Flexibility: many applications and wide range of capacities
- No unpleasant and noxious odours

# **Benefits of Scharlau**

- Outstanding quality
- **♦** Globally available: international sales network
- 30 years experience

# Aquagent® User Guide

	ME0304 Methanol, dry	AQ0011 Aquagent <sup>®</sup> Methanol Fast	AQ0005 Aquagent® Medium K	AQ0009 Aquagent® Buffer	FO0028 Formamide, dry	AQ0002 Aquagent® Solvent	AQ0008 Aquagent® Solvent CM	AQ0010 Aquagent® Solvent OIL	AQ0023 Aquagent <sup>®</sup> Coulometric CG
AQ0007 Aquagent® Complet 2	۵	٥		۵	٥				
AQ0003 Aquagent® Complet 5	٥	٥		٥	٥				
AQ0004 Aquagent® Complet 5K									
AQ0006 Aquagent® Titrant 2									
AQ0001 Aquagent® Titrant 5									
AQ0022 Aquagent® Coulometric A									
AQ0025 Aquagent® Coulometric Oil									

AQ0024 Aquagent® Coulometric AG		



# **Ordering information**

AQUAGENT® PRODU	JCT RANGE			CAPACITY	ART. NO.
Volumetric				500mL	AQ00070500
			Aquagent® Complet 2	1L	AQ00071000
		Reagents		2,5L	AQ00072500
			Aquagent® Complet 5	500mL	AQ00030500
				1L	AQ00031000
				2,5L	AQ00032500
			Aquagent® Complet 5K	500mL	AQ00040500
				1L	AQ00041000
		Solvents	Methanol, dry (max. 0,005% H <sub>2</sub> O), reagent grade	1L	ME03041000
				2,5L	ME03042500
			Aquagent® Methanol Fast	1L	AQ00111000
				2,5L	AQ00112500
				500mL	AQ00050500
			Aquagent® Medium K	1L	AQ00051000
		Additives	Aquagent® Buffer, acid	500mL	AQ00090500
				1L	AQ00091000
			Formamide, dry (max. 0,02% H <sub>2</sub> O), reagent grade	1L	FO00281000
		Titrants	Aquagent® Titrant 2	500mL	AQ00060500
				1L	AQ00061000
			Aquagent® Titrant 5	500mL	AQ00010500
				1L	AQ00011000
				2,5L	AQ00012500
		Solvents	Aquagent® Solvent	1L	AQ00021000
				2,5L	AQ00022500
			Aquagent® Solvent CM	1L	AQ00081000
				2,5L	AQ00082500
			Aquagent® Solvent OIL	1L	AQ00101000
	Cells with diaphragm		Aquagent® Coulometric A, anolyte	500mL	AQ00220500
			Aquagent® Coulometric Oil, anolyte	100mL	AQ00250100
Coulometric			Aquagent® Coulometric CG, catholyte	10 x 5mL	AQ00230050
Coulometric				100mL	AQ00230100
	Cells without diaphragm		Aquagent® Coulometric AG	500mL	AQ00240500
				1L	AQ00241000
	Liquids Solids		Aquagent® standard solution 1.0	10 x 4mL	AQ00190040
Standards			Aquagent® standard solution 10.0	10 x 8mL	AQ00200080
			Aquagent® standard solution 5.0  Aquagent® di-Sodium tartrate dihydrate	100mL	AQ00210100
				500mL	AQ00210500
				25 g	AQ00300025
				100 g	AQ00300100

#### Scharlab S.L.

Gato Pérez, 33. Pol. Ind. Mas d'en Cisa. 08181 Sentmenat, Barcelona, Spain Tel.: +34 93 715 19 40 - Fax: +34 93 715 27 65 E-mail: scharlab@scharlab.com

#### Scharlab Italia S.r.l.

Via Alcide De Gasperi 56. 20070 Riozzo Di Cerro al Lambro (Mi), Italy Tel.: +39 02 9823 0679 / +39 02 9823 6266 Fax: +39 02 9823 0211 / +39 02 9811 9288 E-mail: customerservice@scharlab.it

#### Scharlab Magyarország Kft.

4034 Debrecen, Vágóhíd. u. 2., Hungary Tel: 0036(88)787-634 - Fax: 0036(88)781-081 E-mail: info@scharlab.hu

#### Scharlab Philippines, Inc.

18G Miller Compound, Barangay Bungad, San Francisco Del Monte, Quezon City 1105 Philippines Tel.: +63 2 3514972 - Fax: +63 2 3514972 E-mail: infophilippines@scharlab.ph

#### Scharlab Brasil SA

Estrada do Campo Limpo, 780 São Paulo Tel.: (11) 5512 5744 - Fax: (11) 5511 9366 E-mail: scharlab@scharlab.com.br



Complet Aquagent® product information

