

Get Quality Hand Protection At **New Reduced Prices!**

Ansell's vast range of safety gloves, including Microflex® and Micro-Touch® are designed to provide superior protection from hazardous and infectious substances, delivering the ultimate in hand protection across a variety of applications.

MICRO-FLEX Ergonomic Nitrile Textured Fingertips

Blue Lightweight, Yet strong Gloves

Gentle on sensitive, frail skin due to low residual protein content
Comfortable with high tactile sensitivity



\$28 Box of 250 SAVE 30% on 2022 prices
XL 230 per box

Size	XS (5.5-6)	S (6.5-7)	M (7.5-8)	L (8.5-9)	XL (9.5-10)
Code	ANSE93-833/5.5-6	ANSE93-833/6.5-7	ANSE93-833/7.5-8	ANSE93-833/8.5-9	ANSE93-833/9.5-10

MICRO-FLEX Nitrile Textured Fingertips

Purple Robust Exam Gloves

Textured fingertips ensure better grip when handling tools and equipment
Excellent barrier protection against potentially hazardous substances



\$18 Box of 100 SAVE 45% on 2022 prices

Size	XS (5.5-6)	S (6.5-7)	M (7.5-8)	L (8.5-9)	XL (9.5-10)
Code	ANSE93-843/5.5-6	ANSE93-843/6.5-7	ANSE93-843/7.5-8	ANSE93-843/8.5-9	ANSE93-843/9.5-10

MICRO-FLEX Nitrile + Neoprene, Ambidextrous

Green Chemical-resistant Gloves

Three layers ensure superior protection against harsh chemicals
Soft material and ergonomic design offer outstanding fit, feel and flexibility



\$24 Box of 50 SAVE 34% on 2022 prices

Size	XS (5.5-6)	S (6.5-7)	M (7.5-8)	L (8.5-9)	XL + *XXL
Code	ANSE93-260/5.5-6	ANSE93-260/6.5-7	ANSE93-260/7.5-8	ANSE93-260/8.5-9	ANSE93-260/9.5-10 *260/10.5-11

MICRO-TOUCH® Nitra-Free Nitrile

Pink Medical Exam Gloves

Excellent resistance to a broad range of commonly used chemicals
Protects from latex Type I allergy



\$14 Box 100 SAVE 64% on 2022 prices

Size	S (6.5-7)	M (7.5-8)	L (8.5-9)	XL (9.5-10)
Code	ANSE4511	ANSE4512	ANSE4513	ANSE4514

MICRO-TOUCH SensiClean II

White Latex Multipurpose Gloves

Gentle on sensitive, frail skin due to low residual protein content
Comfortable with high tactile sensitivity



\$14 Box of 100 SAVE 63% on 2022 prices
XL 80 per box

Size	XS (5.5-6)	S (6.5-7)	M (7.5-8)	L (8.5-9)	XL (9.5-10)
Code	ANSE4566	ANSE4560	ANSE4562	ANSE4564	ANSE4568

MICRO-TOUCH DermaClean

White Latex Non-sterile Gloves

Textured outer surface provides an excellent grip in wet and dry conditions
Superior elasticity provides great comfort and flexibility



\$14 Box of 100 SAVE 45% on 2022 prices
XL 80 per box

Size	XS (5.5-6)	S (6.5-7)	M (7.5-8)	L (8.5-9)	XL (9.5-10)
Code	ANSE4570	ANSE4572	ANSE4574	ANSE4576	ANSE4578



Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in NZD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.

Product Information	MICRO-TOUCH SensiClean II	MICRO-TOUCH DermaClean	MICROFLEX 93-833	MICRO-TOUCH NITRAFREE	MICROFLEX 93-843	MICROFLEX 93-260
Product Code	ANSE456X	ANSE457X	ANSE93-833/X	ANSE451X	ANSE93-843/X	ANSE93-2260/X
Colour	White	White	Blue	Pink	Purple	Green
Category	Multipurpose	Multipurpose	High Touch	High Touch	Robust	High Risk
Pack Size - Gloves/Box	100 (XL-80)	100 (XL-80)	250 (XL-230)	100	100	50
Powder Content	Powder Free	Powder Free	Powder Free	Powder Free	Powder Free	Powder Free
“Freedom from Holes (Inspection Level I)”	1.5 AQL	1.5 AQL	.65 AQL	1.5 AQL	.65 AQL	.65AQL
External Glove Surface	Fully Textured	Fully Textured	Textured Fingers	Textured Fingers	Textured Fingers	Textured Fingers
Finger Thickness	0.15mm	0.16mm	0.11mm	0.14mm	0.18mm	0.20mm
Palm Thickness	0.14mm	0.14mm	0.07mm	0.08mm	0.11mm	0.198mm
Cuff Length	Standard	Standard	Standard	Standard	Standard	Extended
Glove Length	240mm	240mm	240mm	240mm	245mm	285mm
Sterile	No	No	No	No	No	No
Sizes Available	XS, S, M, L, XL	XS, S, M, L, XL	XS, S, M, L, XL	S, M, L, XL	XS, S, M, L, XL	XS, S, M, L, XL, 2XL
Tested for use with Chemotherapy Drugs	No	No	No	Yes*	Yes*	Yes*
Allergy Prevention	Low allergenicity latex reduces risk of contact dermatitis	No	Latex (Type 1)	Latex (Type 1) Chemical (Type IV)	Latex (Type 1)	Latex (Type 1)
Key Features	Gentle on sensitive, frail skin due to low residual protein content. Comfortable with high tactile sensitivity.	High grip in wet & dry conditions. Smooth inner surface for easy donning. Good for general applications.	Thin, strong exam glove with Ergoform Technology. Supports musculoskeletal health. Excellent value Nitrile option for broad range of applications.	Chemical accelerator free, minimises the risk of latex allergy Type IV and allergic contact dermatitis.	Durable Nitrile Exam glove with Advanced Barrier Protection. Thickest palm and finger thicknesses across the standard ranges. Utilise the ANSELL GUARDIAN® Chemical Product Selector Tool to check suitability for your application.	Thinnest chemical- resistant synthetic composite disposable glove. Enhanced tactility and dexterity Three-layer design for superior protection against harsh chemicals. Utilise the ANSELL GUARDIAN® Chemical Product Selector Tool to check suitability for your application.
Price/Glove	0.14c (0.17c XL)	0.14c (0.17c XL)	0.11c (0.12c XL)	0.14c	0.18c	0.48c