



Axygen PCR Strips and Tubes

Samples available to trial on request!

- PP Ultra thin and consistent wall thickness for precise thermal transfer
- Even insertion cap for a tight-fitting leakproof seal to prevent evaporation
- Choice of attached or unattached numbered caps
- Available with a domed, flat frosted or a flat clear cap for qPCR
- Available as single tubes or strips of 4, 8 or 12
- RNase-/DNase-free, Human gDNA free, Nonpyrogenic, PCR inhibition free

**Super
Discount!**



Description	Pack Size	Code	Price (excl. GST)	Sale Price (excl. GST)
PCR Tube 0.2ml Flat Cap Assorted Colours	1,000 Items	AXYGPCR-02-A	\$92.20	\$74.20
PCR Tube 0.2ml Flat Cap Clear	1,000 Items	AXYGPCR-02-C	\$83.10	\$59.60
PCR Tube 0.2ml Flat Cap Maximum Recovery Clear	1,000 Items	AXYGPCR-02-L-C	\$120.40	\$99.80
PCR Tube 0.2ml Domed Cap Assorted Colours	1,000 Items	AXYGPCR-02D-A	\$94.00	\$57.00
PCR Tube 0.2ml Domed Cap Clear	1,000 Items	AXYGPCR-02D-C	\$86.00	\$67.90
PCR Tube 0.5ml Flat Cap Assorted Colours	1,000 Items	AXYGPCR-05-A	\$75.60	\$59.40
PCR Tube 0.5ml Flat Cap Clear	1,000 Items	AXYGPCR-05-C	\$75.30	\$60.20
PCR Tube 0.2ml 8-Strip with Individually Attached Flat Cap Clear	126 Items	AXYGPCR-0208-AF-C	\$183.90	\$105.00
PCR Tube 0.2ml 8-Strip with Individually Attached Domed Cap Clear	126 Items	AXYGPCR-0208-AD-C	\$183.90	\$105.00
PCR Tube 0.2ml 8-Strip Tube Assorted Colours Without Cap	125 Items	AXYGPCR-0208-A	\$174.40	\$144.00
PCR Tube 0.2ml 8-Strip Clear Without Cap	125 Items	AXYGPCR-0208-C	\$164.00	\$116.60
PCR Tube 0.2ml 8-Strip With 8-Strip Domed Cap Clear	125 Items	AXYGPCR-0208-CP-C	\$196.80	\$144.50
PCR Cap Domed 8-Strip Assorted Colours	125 Items	AXYGPCR-02CP-A	\$38.50	\$16.00
PCR Cap Domed 8-Strip Cap Clear	125 Items	AXYGPCR-02CP-C	\$33.80	\$25.60
PCR Cap Flat Writing Surface 8-Strip Clear	125 Items	AXYGPCR-02-FCP-C	\$26.20	\$21.20
PCR Cap R-T PCR 8-Strip Ultra-Clear	125 Items	AXYGPCR-2CP-RT-C	\$57.20	\$45.70
PCR Tube 0.1ml 4-Strip Tube with 4-Strip Cap Clear	1,250 Items	AXYGPCR-0104-C_U	\$724.70	\$548.50
PCR Tube 0.1ml Low Profile 8-Strip Tube with R-T PCR Strip Cap Clear	1,250 Items	AXYGPCR-0108-LP-RT-C	\$1,830.50	\$1,393.20
PCR Tube 0.1ml Low Profile 8-Strip Tube with R-T PCR Strip Cap White	1,250 Items	AXYGPCR-0108-LP-RT-W	\$1,830.50	\$1,393.20

If you would like a free 96 well PCR rack (as shown in picture above) to go with your tubes please ask (available while stocks last).

Most items are available ex stock, please enquire for availability.

Contact us for more information
and to place your orders

Bio-Strategy Ltd.

T: 0800 34 24 66

E: sales.nz@bio-strategy.com

www.bio-strategy.com

shop.bio-strategy.com



VALID UNTIL:
30 September 2021

Ref. Code:
AXYG-PCR

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in NZD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.