

Quality Duran Beakers

Low and high form beakers on sale!

Manufactured from Borosilicate glass 3.3 these Duran beakers have very good chemical resistance, virtual inert behaviour, a high maximum usage temperature, high light transparency and minimal thermal expansion.


Featuring fired on highly durable easy-to-read white ceramic graduations and marking area. Printed with Lot numbers for traceability with downloadable certificates.

**Up to
25% OFF!**

Description	Pack Size	Code	Price (excl. GST)	Sale Price (excl. GST)
Low Form Beakers				
Beaker 5ml Duran Borosilicate 22dx30h Spout	Each	SCOT211060701	\$7.70	\$5.80
Beaker 10ml Duran Borosilicate 26dx35h Spout	Each	SCOT211060804	\$8.80	\$6.70
Beaker 25ml Duran Borosilicate 34dx50h Spout	Each	SCOT211061406	\$6.50	\$5.00
Beaker 50ml Duran Borosilicate 42dx60h Spout	Each	SCOT211061706	\$6.30	\$4.80
Beaker 100ml Duran Borosilicate 50dx70h Spout	Each	SCOT211062402	\$6.10	\$4.70
Beaker 150ml Duran Borosilicate 60dx80h Spout	Each	SCOT211062908	\$6.70	\$5.20
Beaker 250ml Duran Borosilicate 70dx95h Spout	Each	SCOT211063604	\$6.70	\$5.20
Beaker 400ml Duran Borosilicate 80dx110h Spout	Each	SCOT211064103	\$7.90	\$6.00
Beaker 600ml Duran Borosilicate 90dx125h Spout	Each	SCOT211064806	\$8.80	\$6.70
Beaker 800ml Duran Borosilicate 100dx135h Spout	Each	SCOT211065305	\$13.50	\$10.30
Beaker 1L Duran Borosilicate 105dx145h Spout	Each	SCOT211065408	\$13.90	\$10.60
Beaker 2L Duran Borosilicate 132dx185h Spout	Each	SCOT211066301	\$25.90	\$19.80
Beaker 3L Duran Borosilicate 152dx210h Spout	Case of 4	SCOT211066807	\$175.70	\$133.90
Beaker 5L Duran Borosilicate 170dx270h Spout	Each	SCOT211067306	\$82.90	\$63.20
Beaker 10L Duran Borosilicate 217dx350h Spout	Each	SCOT211068602	\$398.00	\$303.40
Tall Form Beakers				
Beaker 50ml Duran Borosilicate 38dx70h Spout	Each	SCOT211161704	\$6.50	\$5.00
Beaker 100ml Duran Borosilicate 48dx80h Spout	Each	SCOT211162409	\$6.50	\$5.00
Beaker 150ml Duran Borosilicate 54dx95h Spout	Each	SCOT211162906	\$6.90	\$5.30
Beaker 250ml Duran Borosilicate 60dx120h Spout	Each	SCOT211163602	\$7.30	\$5.60
Beaker 400ml Duran Borosilicate 70dx130h Spout	Each	SCOT211164101	\$7.90	\$6.00
Beaker 600ml Duran Borosilicate 80dx150h Spout	Each	SCOT211164804	\$9.30	\$7.10
Beaker 800ml Duran Borosilicate 90dx175h Spout	Each	SCOT211165303	\$14.60	\$11.30
Beaker 1L Duran Borosilicate 95dx180h Spout	Each	SCOT211165406	\$14.90	\$11.40
Beaker 2L Duran Borosilicate 120dx240h Spout	Each	SCOT211166308	\$27.30	\$20.90
Beaker 3L Duran Borosilicate 135dx280h Spout	Case of 2	SCOT211166805	\$122.50	\$93.40

Please note the 5ml and 10ml low form beakers are supplied without graduations and without a retrace code.

Please note the 10L low form beaker is NON-DIN/ISO size.



Contact us for more information
and to place your orders

Bio-Strategy Ltd.

T: 0800 34 24 66
E: sales.nz@bio-strategy.com
www.bio-strategy.com
shop.bio-strategy.com



VALID UNTIL:
30 September 2021

Ref. Code:
SCOT-BEAKER

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in NZD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.