

Epredia Cytofunnels for Cytospin CytoCentrifuge

Shandon Cytofunnels are a disposable sample chamber with white filter card and cap.

10% discount on:

- EZ Single Cytofunnel with white filter cards
 - Does not require metal Cytoclips
- Single CytoFunnel with white filter cards
 - Requires use of metal Cytoclips
- White filter cards, 2 hole

EPKMA78710020 EZ Single Cytofunnel

Shandon EZ Single Cytofunnel is a disposable sample chamber with white filter card and cap.

- Metal Cytoclips slide clip is not required
- For single use with all Epredia Cytospin cyto centrifuge applications
- Optimises cell yield
- Produces cost effective, thin layer preparations
- Reduces the risk of exposure to pathogenic samples
- Disposable to eliminate time-consuming decontamination and cleaning
- Decreases the possibility of cross contamination
- Use for sample volumes from 0.1-0.5ml, cell deposition area of 6mm in diameter (28mm²)
- Pre-attached white filter card for high absorbency

EPKM1102548 Single CytoFunnel

Shandon Single Cytofunnel is a disposable sample chamber with white filter card and cap.

- Requires the use of stainless steel Cytoclips
- Fit snugly against a specially marked Cytoslide
- Held in place with a stainless steel Cytoclip
- Cytoclips easily unlock to dispose of the sample chamber and to remove the slide
- 12 assemblies fit into the Cytospin head
- Use for sample volumes from 0.1-0.5ml, cell deposition area of 6mm in diameter (28mm²)
- Pre-attached white filter card for high absorbency



Description	Code	Price (excl. GST)	Sale Price (excl. GST)
EZ Single Cytofunnel with White Filter Cards (Shandon), 200 Pcs	EPKMA78710020	\$782	\$703
Single Cytofunnel with White Filter Cards (Shandon), 500 Pcs	EPKM1102548	\$1,138	\$1,024
White Filter Cards, 2 Hole, 200pcs	EPKM5991022	\$114	\$102

EPKM5991022 White Filter Cards, 2 hole

- Highly absorbent disposable filter cards for use with TPX chamber
- Not for use with Cytofunnel disposable sample chambers
- Each filter card has two holes to ensure proper orientation with the TPX sample chamber
- White filter cards are used for high absorbency



Contact us for more information and to place your orders

Bio-Strategy Ltd.

T: 0800 34 24 66
E: sales.nz@bio-strategy.com
www.bio-strategy.com
shop.bio-strategy.com

 [linkedin.com/company/bio-strategy](https://www.linkedin.com/company/bio-strategy)

VALID UNTIL:
30 September 2021

Ref. Code:
EPRN-CONS

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in NZD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.