

Epredia Cytospin 4 Cyto centrifuge

***Order a CytoSpin 4 and receive 50% off a starter kit.**

The Gold standard in Cyto centrifugation

- Reliable
- Easy to use
- Multiple applications: Cytology, Microbiology, Haematology/ Oncology, and Research
- Use with Epredia consumables for best results

The Cytospin™ 4 cyto centrifuge has been synonymous with cyto centrifugation for over 40 years and is used in laboratories around the world.

When combined with Epredia cytology consumables, it provides economical thin-layer preparations from any liquid matrix, especially hypocellular fluids such as spinal fluid and urine.

Convenience and speed spun into one

- Process twelve specimens at one time, allowing for several protocols at once
- Accepts all protocols from previous generations of Cytospin cyto centrifuges
- Lid release mechanism allows for one-handed opening and closing
- Removable program card allows for up to nine programs to be logged for fast and convenient reference, and protected program memory stores up to 23 routines for instant recall
- Control panel allows user to control power and time, and monitor speed
- Specimen safety alarm feature reminds users at one-minute intervals to remove specimens, protecting them from air-drying

Designed for user safety

- Only operates when the lid is locked and remains locked during rotation
- Polycarbonate window allows the operator to see the sealed head during cyto centrifugation



**50% Off
Starter Kit***

Description	Code	Price (excl. GST)	Sale Price (excl. GST)
Cytospin 4 Unit with Rotor	EPRNA78300002	\$18,802	
Clipped Funnel Starter Kit (with Reusable Metal Clips)	EPRN3120110	\$2,250	\$1,125
EZ Funnel Starter Kit (all Plastic/Disposable)	EPKMA78710013	\$1,062	\$531

Simple Principle for Excellent Presentation

- Using centrifugal principles, the cyto centrifuge deposits cells onto a clearly defined area of a glass slide and allows for the absorption of the residual fluid into the sample chamber's filter card
- Cyto centrifugation constructively flattens cells for excellent nuclear presentation
- The instrument's spinning action tilts the Cytofunnel sample chamber upright and centrifuges cells onto the deposition area of the slide, creating a cellular monolayer that provides all cell types equal opportunity for presentation
- In the "load" or "stop" position, this unique tilting feature reduces cell loss by preventing residual fluid from coming into contact with the prepared slide

Starter Kits for Cytospin 4 Cyto centrifuge

EPRN3120110 Clipped funnel Starter Kit

Kit with reusable metal clips

Includes:

- 50x Single White Cytofunnels
- 25x Double Cytofunnels
- Cytoslides
- 12 Cytoclips (reusable metal clips)

EPKMA78710013 EZ funnel Starter Kit

Disposable Plastic Kit

Includes:

- 25x EZ Megafunnels
- 20x EZ Single White Cytofunnels
- 20x EZ Double Cytofunnels
- Cytoslides

POWERED BY SHANDON MICROM MENZEL-GLASER Richard Allan Scientific

Contact us for more information
and to place your orders

Bio-Strategy Ltd.

T: 0800 34 24 66
E: sales.nz@bio-strategy.com
www.bio-strategy.com
shop.bio-strategy.com

 [linkedin.com/company/bio-strategy](https://www.linkedin.com/company/bio-strategy)

VALID UNTIL:
30 September 2021

Ref. Code:
EPRN-CYTO4

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in NZD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.