

## Leica DM750 Microscope

### Order a DM750 and receive 15% discount off ICC50W camera.

#### The Leica DM750 Microscope for high quality images and reliable performance!

Leica Optics + Bright LED illumination leads to a higher quality image.

- Optics are based on the same optical platform as Leicas' research microscopes
- EZLite™ – provides over 20 years of bright, LED illumination which saves time and the cost of replacement bulbs
- Auto off – patented time delay shutoff saves energy

Robust: Quality materials and construction gives reliable performance over a longer life.

- Weighted focus knobs with internal brass mechanism for accurate focus capability
- Wear resistant, special stage finish offers additional protection from friction damage
- AGTREAT™ on all microscope touchpoints, inhibits bacteria growth and helps prevent the spread of disease from person to person
- EZStore™ design with handle and cord wrap allows easy carrying, easy lifting and easy cord storage

#### Would you like to share, capture and archive your microscope images?

Add: ICC50W High-Definition Wireless Camera module.

5 modes/choices of connectivity and working with images:

- WiFi mode + Leica AirLab App for mobile devices
- USB mode connect to PC and Leica software
- HDMI port direct to HD screens
- SD card
- Ethernet

#### DM750 microscope – customise to your needs and budget, starts from \$2,300.



#### Special offer: Leica DM750 Microscope ex-demo (1x only)

Brightfield and phase contrast (Phase slider).

Koehler illumination. Objectives: HI Plan 4x, 10x, 20x, 40x, and N plan 100x oil. 10x focussing eyepieces.

Description	Code	Price (excl. GST)	Sale Price (excl. GST)
Ex-demo: DM750 Microscope with Phase Contrast (1x only)	LEIM-DM750	\$8,412	\$5,468
Leica ICC50W camera Module	LEIM-13613735	\$2,372	\$2,016
DM750 Microscope	LEIM-XXX	<b>Starts from \$2,300</b>	

Contact us for more information and to place your orders

#### Bio-Strategy Ltd.

T: 0800 34 24 66  
E: sales.nz@bio-strategy.com  
www.bio-strategy.com  
shop.bio-strategy.com

 [linkedin.com/company/bio-strategy](https://www.linkedin.com/company/bio-strategy)

VALID UNTIL:  
30 September 2021

Ref. Code:  
LEIM-DM750

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in NZD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.