

Miele PG8583 CD Washer/Disinfector

10% package discount including free detergent, neutraliser and commissioning!

This underbench /standalone model combines a powerful pump with integrated heating elements for safe washing cycles up to 93°C. The DryPlus active HEPA filtered drying ensures loads are dry and ready to use.

- 16 programs including new "Mini" program for lightly soiled glassware
- Conductivity monitoring – proof that load is "clean"
- Configured for washing wide neck and narrow neck glassware
- Dispensing pumps for liquid detergent and liquid neutraliser
- Compact: 900mm x 700mm x 835mm W x D x H (820mm H with lid removed)
- Pullout drawer for housing chemical canisters (detergent and neutraliser)
- Stainless steel interior/exterior for durability and cleanliness
- ADP model including pump for AD water and waterproofing kit
- Many other options and accessories to select from
- Spray arm and spray pressure monitoring

Package Deal!



Description	Code	Price (excl. GST)	Sale Price (excl. GST)
Miele PG8583 CD Package 1 Includes: A102 Upper Basket with 2 x AK12 Inserts, A150 Lower Module Basket with 2 x A301/5 Modules (3 x 6 = 18 Nozzles/Module), Waterproofing Kit and commissioning	MIELPG8583CD	\$29,950	\$26,950

Contact us for more information
and to place your orders

Bio-Strategy Ltd.

T: 0800 34 24 66
E: sales.nz@bio-strategy.com
www.bio-strategy.com
shop.bio-strategy.com



[linkedin.com/company/bio-strategy](https://www.linkedin.com/company/bio-strategy)

VALID UNTIL:
30 September 2021

Ref. Code:
MIEL-PG8583CD

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in NZD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.