

## Priorclave Autoclave Range

**Purchase a Priorclave autoclave and receive a free Priorclave thermal printer valued at over \$2,000.**

**British manufacturer Priorclave have an autoclave model tailored to suit your requirements.**

- Cylindrical chamber front loading models with 40L, 60L, 100L, 150L, 200L, 320L, and 400L
- Cylindrical chamber top loading models with 60L, 85L, 100L, and 150L
- Rectangular section chamber models with 230L, 285L, 350L, 450L and, 700L
- Manufactured to stringent British standards the autoclaves are compliant to AS/NZS standards and thermal interlocks ensure operator safety preventing access to the chamber until load is cooler than 80°C
- Wide range of accessories including vacuum options, air ballasting, water jackets, water baffles, wandering probe, multiple programs, trolleys, loading cassettes, discard containers, printer, water softener, drain condenser, steam generators, filters, hoist etc.
- Manual or power door as well as pass through models
- Grade 316SS chamber with 10-year warranty
- Top loading units supplied with two baskets and front loading units with two shelves as standard
- All models supplied complete with USB port, timed free-steaming, fan assisted cooling and media warming facility

**Contact us now with your autoclave requirements.**

Description	Code	Price (excl. GST)	Sale Price (excl. GST)
Priorclave Thermal Printer with Print Delay Feature	PRIOPC/PRN/000	\$2,250	<b>Free</b>



**Free  
Printer!**



Contact us for more information  
and to place your orders

**Bio-Strategy Ltd.**

T: 0800 34 24 66

E: sales.nz@bio-strategy.com

www.bio-strategy.com

shop.bio-strategy.com

 [linkedin.com/company/bio-strategy](https://www.linkedin.com/company/bio-strategy)

VALID UNTIL:  
30 September 2021

Ref. Code:  
PIRO-AUTO

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in NZD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.