

PHCbi Compact -80°C Chest Freezer, 84 litres

With Covid 19 impacting on our population – freezers for vaccine storage are in high demand.

Ideal low temperature compact freezer for vaccine storage.

- 84 litres with small footprint: 550 x 685mm
- Adjustable temperature range from -60°C to -80°C
- Single compressor system, reliable with three year warranty
- Excellent brand: PHCbi (formerly Panasonic / Sanyo)
- Filter-less design: reduced maintenance
- Alarms to protect valuable product
- Low noise design: 47dB(A)
- VIP Plus vacuum insulation panel technology
- Improved energy consumption achieving energy savings of approx. 40%
- Dry remoter alarm contacts and access port
- Suitable for operation on 240V 50Hz 10 amp supply



**20%
Off!**

Description	Code	Price (excl. GST)	Sale Price (excl. GST)
PHCbi -80°C Compact Chest Freezer, 84 litres	SANYMDF-C8V1-PE	\$9,517	\$7,614

Contact us for more information
and to place your orders

Bio-Strategy Ltd.

T: 0800 34 24 66
E: sales.nz@bio-strategy.com
www.bio-strategy.com
shop.bio-strategy.com

 [linkedin.com/company/bio-strategy](https://www.linkedin.com/company/bio-strategy)

VALID UNTIL:
30 September 2021

Ref. Code:
SANY-84L

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in NZD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.