

Axygen® PCR Microplates

Stock up for the Holiday Period!

Axygen 96 well and 384 well clear and white PCR microplates with black alpha-numerical lettering allow for easy observation and tracking. Axxygen PCR microplates are manufactured to stringent specifications under strict process controls and are assembled via automation to ensure consistent performance and quality. The microplates have a thin and consistent wall thickness for precise thermal transfer along the entire surface of the microplate wells.


Features and Benefits

- 96 well and 384 well
- Comes in clear and white for standard and Real-Time PCR or R-T PCR
- Available in a wide variety of sizes and configurations
- Compatible with most of the major thermal cyclers
- Automation friendly
- Raised well rims
- Manufactured from a homopolymer with USP Class 6 virgin polypropylene
- Certified RNase-/DNase-free
- Nonpyrogenic

Applications

- Standard, gradient, and real-time PCR
- High throughput PCR

**STOCK UP
& SAVE!**



Code	Description	Pack Size	List Price	Sale Price
AXYGPCR-24-C	Plate PCR 24 well PP clear	50 items	\$133.60	\$77
AXYGPCR-32-C	Plate PCR 32 well PP clear	50 items	\$183.10	\$121
AXYGPCR-96-AB-C	Plate PCR 96 well clear semi skirted (#1; #7)	50 items	\$404.90	\$301
AXYGPCR-96M2-HS-C	Plate PCR 96 well clear half skirt with single notch (#1; #3; #4; #6; #7; #9; #10)	50 items	\$316.90	\$249
AXYGPCR-96-LP-AB-C	Plate PCR 96 well clear half skirt (#2; #8)	100 items	\$345.30	\$237
AXYGPCR-96-LP-FLT-C	Plate PCR 96 well clear half skirt low profile for ABI (#3; #10)	100 items	\$439.80	\$346
AXYGPCR-96-LP-FLT-W	Plate PCR 96 well white half skirt low profile for ABI (#3)	100 items	\$589.10	\$438
AXYGPCR-96-FS-C	Plate PCR 96 well clear full skirt (#3; #4; #9; #10)	50 items	\$321.70	\$240
AXYGPCR-96-FLT-C	Plate PCR 96 well clear no skirt (#4; #6; #9; #10)	100 items	\$507.40	\$349
AXYGPCR-96-LC480-C-NF	Plate PCR 96 well clear for Roche 480 (#5)	50 items	\$474.80	\$353
AXYGPCR-96-LC480-W-NF	Plate PCR 96 well white PP for Roche 480 (#5)	50 items	\$281.50	\$210
AXYGPCR-96-HS-C	Plate PCR 96 well clear half skirt thin wall double notch corner	50 items	\$369.40	\$293
AXYGPCR-96-C	Plate PCR 96 well clear (#6)	50 items	\$323.00	\$249
AXYGPCR-96-SG-C	Plate PCR 96 well clear segmented and pre scored for Beckman CEQ	50 items	\$315.40	\$235
AXYGPCR-384M2-C	Plate PCR 384 well clear skirted for ABI (#1; #7; #9; #10)	50 items	\$273.20	\$204
AXYGPCR-384-LC480-W-NF	Plate PCR 384 well white for Roche 480 (#5);	50 items	\$438.90	\$327
AXYGPCR-384-LC480WNFBC	Plate PCR 384 well white with bar codes for Roche 480 (#5);	50 items	\$608.70	\$479

More products on next page

Axygen® PCR Microplates

Code	Description	Pack Size	List Price	Sale Price
Barcoded Plates				
AXYGPCR-96-LP-AB-C-PBC	Plate PCR 96 well clear half skirt with bar codes (#2; #8)	100 items	\$535.10	\$421
AXYGPCR-96-AB-C-PBC	Plate PCR 96 well clear semi skirted with bar codes (#1; #7)	50 items	\$592.50	\$466
AXYGPCR-96-FS-C-PBC	Plate PCR 96 well clear full skirt with bar codes (#3; #4; #9; #10)	50 items	\$385.60	\$303
AXYGPCR-96-LC480-W-BC	Plate PCR 96 well white PP with bar codes and sealing film for Roche 480 (#5)	50 items	\$929.60	\$731
AXYGPCR-384M2-C-PBC	Plate PCR 384 well clear skirted with bar codes for ABI (#1; #7; #9; #10)	50 items	\$526.20	\$414
AXYGPCR-384-LC480WNFBC	Plate PCR 384 well white with bar codes for Roche 480 (#5)	50 items	\$608.70	\$479

NB: Due to Covid-19 some PCR bar coded plates are experiencing backorder delays, please inquire for availability.

Axygen Axyspin Plate Spinner




The Axigen Axyspin Plate Spinner is the first centrifuge designed specifically for quick and easy centrifugation of samples in PCR plates or microplates.

- Capacity for two standard plates
- Top loading design
- Exceptionally small footprint
- Starts and stops in seconds



Code	Description	List Price	Sale Price
AXYGPLATESPINNER-230AU	Centrifuge PCR miniplate spinner Axigen 2 plate 230v	\$1,244.90	\$979

Axygen Sealing Film



Code	Description	Pack Size	List Price	Sale Price
AXYGUC-500	Sealing film ultraclear optical RT-PCR polyolefin 141x77 (#1; #2; #3; #4; #6)	100 items	\$275.40	\$205
AXYGPCR-TS	Sealing film cyclerseal PCR/Storage polyester 135x80mm (#8; #9; #10)	100 items	\$102.50	\$79
AXYGPCR-AS-200	Sealing film aluminium pierceable PCR/Storage 146x79mm	100 items	\$204.30	\$161
AXYGPCR-SP	Sealing film Axysel ELISA/Tissue Culture polyester 146x79mm	100 items	\$80.00	\$59

Contact us for more information and to place your orders

Bio-Strategy Ltd.

T: 0800 34 24 66
E: sales.nz@bio-strategy.com
www.bio-strategy.com
shop.bio-strategy.com



[linkedin.com/company/bio-strategy](https://www.linkedin.com/company/bio-strategy)

VALID UNTIL:
31 December 2021

Ref. Code:
AXYG-PCR

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in NZD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.