

## Epredia PrintMate AS Cassette Printer

---

### Clean Print, Lean Process

**Increase efficiency while reducing the likelihood of errors in your labeling workflow with the PrintMate AS™ cassette printer.**

- Save time – prints each cassette in ten seconds or less
- Compact footprint 38cm x 43.5cm (DxW)
- Print only the cassettes you need, right at the point of use
- Eliminate errors endemic to handwriting processes

### Enhance Efficiency

Achieving a Lean workflow requires having the tools you need where you need them. The PrintMate AS cassette printer's footprint fits easily next to the grossing area. Its printing speed allows grossing at your pace, without long pauses waiting for a labeled cassette.

### Increase Accuracy

The PrintMate AS printer's high-quality printing system marks every cassette with easy-to-read, permanent black type that enhances your lab's accuracy and sample integrity, allowing you to focus on what matters most - positive patient outcomes.

### Adapt to Your Needs

The PrintMate AS cassette printer can be used on-demand in any size laboratory. An optional dual-track collection system collects up to 100 cassettes in organizer trays which stack easily and simplify matching cassettes to specimen containers.



Description	Code	Price
PrintMate AS 450 (6 hopper)	EPRNA8380002	from \$35,700

*More information on next page*

## Epredia Tissue Cassettes for PrintMate Printers

**Epredia tissue cassettes are engineered to optimize print quality when used with PrintMate AS cassette printers.**

- Redesigned for better performance
- Comprehensively validated for print quality, processing quality and reliability
- Tubes are precisely engineered to reduce jamming and prevent downtime

**Buy More  
&  
Save**

**Ask us about standing order pricing!**

Description	Code	Price
Embedding Cassette, round pore style, Lilac, no lids, 750pcs	EPKMB1000732LI	\$181
Embedding Cassette, round pore style, White, no lids, 750pcs	EPKMB1000732WH	\$150
Cassette 2, Slotted, White, with Loose-Packed Lids, 750pcs	EPKMB851729WH	\$141
Histoscreen Cassette, pink, loose with loose packed lids 1000PCS	EPKMB851000PK	\$776



Embedding Cassette  
round pore style

Histoscreen Cassette  
mesh pore style



Cassette 2  
slotted pore style

*More cassette types and colours available*

The entire range of cassettes from Epredia™ has been redesigned for better performance – from printing through to laboratory processes.

### Printing

- The all-new surface finish and reinforced print faces are engineered to optimize print quality when used with PrintMate™ AS cassette printers from Epredia
- All-new PrintMate AS-compatible tubes are precisely engineered to reduce jamming and prevent downtime

### Grossing, Embedding and Processing

- Opening and closing forces have been carefully adjusted for ease of use, while safeguarding valuable specimens during processing
- Design and flow characteristics tested to ensure efficient processing

### Sectioning

- The new cassettes are precision engineered for a solid connection to Epredia microtomes

### Archiving

- The new surface texture is designed to optimize scanability when used with the Arcos™ tissue block management system from Epredia

### Epredia can offer a range of validated cassettes

- All cassette designs have been comprehensively validated for print quality, processing quality and reliability

Contact us for more information  
and to place your orders

**Bio-Strategy Ltd.**

T: 0800 34 24 66

E: [sales.nz@bio-strategy.com](mailto:sales.nz@bio-strategy.com)

[www.bio-strategy.com](http://www.bio-strategy.com)

[shop.bio-strategy.com](http://shop.bio-strategy.com)



**bio-strategy**  
delivering technology

VALID UNTIL:  
31 December 2021

Ref. Code:  
EPRN-MATE



[linkedin.com/company/bio-strategy](https://www.linkedin.com/company/bio-strategy)

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in NZD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.