

## ESCO NordicSafe Class II Cabinet

**Free adjustable height stand, UV light and 2 power sockets (valued at >\$1,800)**

**One of the quietest cabinets on the market today!**

Providing operator, product and environment protection, the NC2 models comply to the EN12469 standard.

The motorised sash automatically activates internal lighting and airflow.

- Quickstart mode- turns on cabinet, light and fan with single button press
- Interlocked UV light with timer – no accidental burns
- 3-year warranty
- Ergonomic sloping front combined with tempered glass side panels allows entry of more natural light
- ULPA downflow and exhaust filters with 99.999% efficiency, 10X better than HEPA filter (99.99%) improves safety and delivers ISO Class 3 air cleanliness
- Sentinel Gold microprocessor controller with digital display shows cabinet status
- Multi-piece stainless steel work surface is easy to access and clean
- Half Speed Mode reduces energy consumption by a further 60%
- Integrated temperature compensated airflow monitoring system
- Extremely low energy consumption: 187 / 272 watts respectively
- Excellent white light illumination: 1,200 / 1,600 watts respectively

**15% Off  
+  
Free Extras!**



Description	Code	List	Price
1.2m Class II Microbiological Safety Cabinet	ESCONC2-4L8	\$16,290	<b>\$13,845</b>
1.8m Class II Microbiological Safety Cabinet	ESCONC2-6L8	\$19,549	<b>\$16,616</b>

Contact us for more information and to place your orders

**Bio-Strategy Ltd.**

T: 0800 34 24 66

E: sales.nz@bio-strategy.com

www.bio-strategy.com

shop.bio-strategy.com



**bio-strategy**  
delivering technology



[linkedin.com/company/bio-strategy](https://www.linkedin.com/company/bio-strategy)

VALID UNTIL:  
31 December 2021

Ref. Code:  
ESCO-NC2

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in NZD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.