

CYTIVA Vivaspin™ Protein Concentrator Spin Columns

Top quality products on sale, order today!

Vivaspin™ ultrafiltration spin columns are designed for fast, non-denaturing concentration and desalting of proteins and other biological samples.

- High recovery rates – up to 30-fold concentration of the sample can be achieved, with recovery of the target molecule typically exceeding 95%
- Sample protection – the patented Vivaspin dead-stop technology ensures samples cannot be concentrated to dryness and enables direct concentrate recovery
- Ease-of-use – one-step sample concentration in a single tube, for minimal sample handling and reduced sample loss
- pH compatibility – for a broad range of pH 1 to 9
- High flow rates – and minimal membrane blockage through the vertical polyethersulfone (PES) membrane design
- Minimal protein adsorption – thanks to the Vivaspin polycarbonate tube design

10% OFF



Description	Code	List	Sale
CYTIVA CONCENTRATOR MWCO			
3000 VIVASPIN 500 [25 items]	GEHE28-9322-18		
5000 VIVASPIN 500 [25 items]	GEHE28-9322-23		
10000 VIVASPIN 500 [25 items]	GEHE28-9322-25	\$148	\$133.20
30000 VIVASPIN 500 [25 items]	GEHE28-9322-35		
100000 VIVASPIN 500 [25 items]	GEHE28-9322-37		
CYTIVA CONCENTRATOR MWCO			
3000 VIVASPIN 6 [25 items]	GEHE28-9322-93		
10000 VIVASPIN 6 [25 items]	GEHE28-9322-96		
30000 VIVASPIN 6 [25 items]	GEHE28-9323-17	\$309	\$278.10
50000 VIVASPIN 6 [25 items]	GEHE28-9323-18		
100000 VIVASPIN 6 [25 items]	GEHE28-9323-19		

Description	Code	List	Sale
CYTIVA CONCENTRATOR MWCO			
3000 VIVASPIN 20 [12 items]	GEHE28-9323-58		
5000 VIVASPIN 20 [12 items]	GEHE28-9323-59		
10000 VIVASPIN 20 [12 items]	GEHE28-9323-60		
30000 VIVASPIN 20 [12 items]	GEHE28-9323-61	\$219	\$197.10
50000 VIVASPIN 20 [12 items]	GEHE28-9323-62		
100000 VIVASPIN 20 [12 items]	GEHE28-9323-63		
CYTIVA CONCENTRATOR MWCO			
3000 VIVASPIN 2 [25 items]	GEHE28-9322-47		
10000 VIVASPIN 2 [25 items]	GEHE28-9322-58	\$235	\$211.50
100000 VIVASPIN 2 [25 items]	GEHE28-9322-58		

Contact us for more information and to place your orders

Bio-Strategy Ltd.

T: 0800 34 24 66
E: sales.nz@bio-strategy.com
www.bio-strategy.com
shop.bio-strategy.com

 [linkedin.com/company/bio-strategy](https://www.linkedin.com/company/bio-strategy)

VALID UNTIL:
31 December 2021

Ref. Code:
GEHE-VIVA

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in NZD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.