

Gyrozen 416 Centrifuge

Compact Low Speed Centrifuge

Compact & versatile, the Gyrozen 416 Centrifuge can easily be configured to meet your requirements!

Commonly used in the fractionation of blood samples and simple precipitation of biological samples.

A separate button for "Soft" start/stop function to achieve very gentle acceleration and deceleration.

Furnished with unique air flow design and anti-vibration structure to minimize heat and noise generation.

Lid drop protection design allows easiness and safety when loading or unloading samples.

Automatic lid release and vivid end alarm for easy retrieval of sample.

A large assortment of rotors, buckets and adaptors for microtubes and all sizes of tubes up to 100 ml.

- Max. Speed: 4,000rpm ; Max Force: 2,774 x g
- 10 program memory & selectable "Soft Start"
- Many options / accessories incl. fixed angle/ swing-out rotors + microplate carriers
- Key safety features: imbalance cut out and lid lock protection
- Maximized capacity of 400ml
- 2-year warranty



Package Deal!

Description	Code	List	Sale
Gyrozen Compact Low Speed Centrifuge Package	GYRO-GZ-0416-PK	\$5,262	\$4,605

Package Deal Includes:

- 1 x 416 Compact Low Speed Centrifuge
- 1 x Rotor Swing Out 4x100ml
- 1 x 100ml Bucket + Biosafety Cap
- 1 x Adaptors*

*Contact Bio-Strategy for adaptor options

Contact us for more information and to place your orders

Bio-Strategy Ltd.

T: 0800 34 24 66
E: sales.nz@bio-strategy.com
www.bio-strategy.com
shop.bio-strategy.com