

Gyrozen 1312 Mini Microcentrifuge

- 12 microtubes or 10 cryotubes can be run up to 13,500 rpm (12,300 xg)
- Compact design with patented air-cooling to minimise noise and heat generation
- Supplied with autoclavable 12-place fixed angle micro-rotor, snap-fit plastic rotor lid and 12 adaptors for 0.2 ml PCR tubes
- A separate "PULSE" key for quick spin



Gyrozen 406 Low Speed Centrifuges

- Compact and robust small centrifuge 406 is ideal for routine low-speed centrifugation
- Commonly used in the fractionation of blood samples and simple precipitation of biological samples
- Automatic lid release and vivid end alarm for easy retrieval of samples
- Supplied with an angle rotor including 6 of aluminum sleeves for 15 ml tubes



Gyrozen 416 Low Speed Centrifuges

- Compact bench top centrifuge with maximized capacity of 400 ml
- 10 program memory & selectable "Soft Start"
- Furnished with unique air flow design and anti-vibration structure to minimize heat and noise generation
- Lid drop protection for convenient and safe sample loading/unloading
- Automatic lid release and vivid end alarm for easy retrieval of samples
- A large assortment of rotors, buckets and adaptors for microtubes and all sizes of tubes up to 100 ml



Gyrozen 1536 Micro Centrifuges

- Compact design built to accommodate PCR tubes, microtubes and cryotubes (0.2 to 2.0 ml)
- Trendy touch wheel and soft touch control panel with key lock function
- 3.5" LCD window with simultaneous display of set and actual parameters, which can be changed during operation
- A separate "PULSE" key for quick spin
- "Soft touch" motorized lid mechanics, with built-in lid drop protection



Description	Quantity	Code	Price NZ\$	Sale NZ\$
1312 Mini Microcentrifuge, white with fixed angle rotors and 12 x 0.2 ml adaptors	5	GYROGZ-1312	1,211	995
406 Compact Low Speed Centrifuge with 6 x 15 ml fixed angle rotor	2	GYROGZ-0406	1,979	1,610
416 Compact Low Speed Centrifuge with swing-out rotor, 100ml buckets & adaptors	4	GYROGZ-0416-PK	5,262	4,605
1536 Microcentrifuge without rotor (various rotors and adaptors available)	1	GYROGZ-1536	2,777	2,250

Contact us for more information
and to place your orders

Bio-Strategy Ltd.

T: 0800 34 24 66
E: sales.nz@bio-strategy.com
www.bio-strategy.com
shop.bio-strategy.com



bio-strategy
delivering technology



[linkedin.com/company/bio-strategy](https://www.linkedin.com/company/bio-strategy)

Valid until:
31 Dec 2021
or while stocks last

Ref. Code:
GYRO-StockSale

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in NZD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.