

Explorer Balances

- 4 programmable touchless sensors
- 14.5 cm (diagonal) Full-Colour touch screen
- Modular design (detachable screen)
- Automatic calibration
- Up to 14 built-in applications including pipette adjustment, ingredient costing, SQC
- Lockout features, password protection
- 5-year warranty



Description	QTY	Code	Sale NZ\$
EX224/AD - 220G X 0.1MG	1	OHAU30061977	4,995
EX4202 - 4200G X 0.01G	1	OHAU83021357	3,995
EX35001 - 35000G X 0.1G	1	OHAU30057065	5,750

Scout Portable Balances

- Durable and reliable, with superior overload protection
- Fast weighing speed with high resolution ensures excellent reliability and repeatable results
- Practical, stackable design for convenient storage
- Multiple connectivity options
- 1-year warranty



Description	QTY	Code	Sale NZ\$
Scout® STX Touch Screen			
STX223 - 220G X 0.001G	1	OHAU30268893	990
Scout® SKX Education			
SKX222 - 220G X 0.01G	2	OHAU30268917	495
SKX422 - 420G X 0.01G	1	OHAU30268918	665
SKX622 - 620G X 0.01G	1	OHAU30268919	715
SKX621 - 620G X 0.1G	1	OHAU30268923	495

Adventurer Balances

- 10.9 cm diagonal full-colour touch screen
- Automatic calibration
- 9 built-in applications
- ≤ 1.5 sec stabilization time
- Die-cast metal base
- 3-year warranty



Description	QTY	Code	Sale NZ\$
AX324 - 320G X 0.1MG	1	OHAU30100604	2,995
AX223 - 220G X 0.001G	1	OHAU30100605	1,995
AX622 - 620G X 0.01G	1	OHAU30122632	1,715

MB120 Moisture Balances

- ABS top housing, stainless steel pan support, stainless steel metal pan handler, in-use cover
- 7 shut-off criteria (manual, timed, autofree %/s, autofree mg/s, A30, A60 and A90)
- 4 heating profiles (standard, fast, ramp and step)
- 40-230°C heating temperature range (1°C increments)
- 13 operational languages (display only, no key in)
- Storage space for up to 100 methods and 1000 total results
- Three-level user management
- 3-year warranty



Including software SmartGuide - assists in identifying the optimal drying temperature for a given sample; SmartGuide analyses your sample and provides a program recommendation in 30 min

Description	QTY	Code	Sale NZ\$
MB120 MOISTURE ANALYSER	2	OHAU30246943	4,250

Navigator Multi-Purpose Balances

- Large, backlit LCD displays for any working environment
- Produces accurate results within one second
- User-friendly with easy set-up
- With ABS housing, level indicator and stainless-steel pan/platform
- 1-year warranty



Description	QTY	Code	Sale NZ\$
NV2201 - 2200 G X 0.1G	2	OHAU30467933	395
NV422 - 422G x 0.01G	1	OHAU30467928	415

Incubating Rocking Shaker

- Combines smooth rocking motion and general-purpose incubation in one compact bench unit
- Electric tilt adjustment from 0 to 15° while unit is operating
- LED displays for temperature, time, and speed and tilt angle
- Timer with audible alarm
- 1-year warranty



Description	QTY	Code	Sale NZ\$
INCUBATING SHAKER ROCKING - ISRK04HDGAU	1	OHAU30391976	2,950



Aquasearcher Bench pH Meter Packages

- A highly reliable and user-friendly range of bench meters
- Designed to support a variety of laboratory requirements – from advanced research to standard application
- Warranty: 1 year on meters | 6 months on electrodes
- Packages include:
pH meter, pH electrode, electrode stand, 4 x 50ml mini buffer kit



Description	QTY	Code	Sale NZ\$
AQUASEARCHER - AB33PH-F	2	OHAU30589827	1,236

Description	QTY	Code	Sale NZ\$
AQUASEARCHER - AB41PH-F incl. STTEMP probe + AS20 stirrer	1	OHAU30589831	1,724

Starter Pens

- Simple, fast, straightforward, no-frills yet highly accurate operation time and time again
- With durable plastic housing, a protective cap that safeguards the sensor and an automatic shutdown feature to preserve battery life
- 1-year warranty



Description	QTY	Code	Sale NZ\$
ST10 - 0 – 14PH, ±0.1	1	OHAU30073970	99
ST20 - 0 – 14PH, 0.0 – 99.0 °C, ±0.05 pH / ± 0.5 °C	1	OHAU30073971	170

Mini Vortex Mixers - Analogue or Digital Models

- Analogue offering variable speed control, and digital allowing the input of exact speed and time for repeatable results
- User-friendly, robust and durable design
- 2-year warranty



Description	QTY	Code	Sale NZ\$
MINI ANALOG - VXMNALAU	1	OHAU30392118	395
MINI DIGITAL - VXMNDGAU	1	OHAU30392125	575

SF40A Printer

- Easy-to-use and set up dot-matrix printer
- Exceptional value and functionality
- Ideal for use in conjunction with OHAUS balances and scale
- 1-year warranty



Description	QTY	Code	Sale NZ\$
SF40A IMPACT PRINTER	2	OHAU30045641	599

Contact us for more information
and to place your orders

Bio-Strategy Ltd.

T: 0800 34 24 66
E: sales.nz@bio-strategy.com
www.bio-strategy.com
shop.bio-strategy.com



bio-strategy
delivering technology

 [linkedin.com/company/bio-strategy](https://www.linkedin.com/company/bio-strategy)

Valid until:
31 Dec 2021

Ref. Code:
OHAU-StockSale

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in NZD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.