

LIMITED STOCK - ORDER NOW!

Ultra-Low Temperature Freezers



MDF-C8V1-PE



MDF-U55V-PE



MDF-U74V-PE

Description	Quantity	Code	Price NZ\$	Sale NZ\$
MDF-C8V1-PE - Chest Freezer, 84L -80°C, VIP Plus vacuum insulation, single compressor system with approx. 40% reduction in power consumption, hot gas bypass around door seal to reduce ice build-up, filter-less design: no messy cleaning dirty filters!	7	SANYMDF-C8V1-PE	9,517	6,975
MDF-U55V-PE - Upright Freezer, 519L -86°C, VIP Plus vacuum insulation, cascade refrigeration system, hot gas bypass around door seal to reduce ice build-up, insulated inner doors, supplied with 3 adjustable shelves	1	SANYMDF-U55V-PE	18,233	13,500
MDF-U74V-PE - Upright Freezer, 728L -86°C, VIP Plus vacuum insulation, cascade refrigeration system, hot gas bypass around door seal to reduce ice build-up, insulated inner doors, supplied with 3 adjustable shelves	2	SANYMDF-U74V-PE	22,176	15,995

Contact us for more information
and to place your orders

Bio-Strategy Ltd.

T: 0800 34 24 66
E: sales.nz@bio-strategy.com
www.bio-strategy.com
shop.bio-strategy.com



bio-strategy
delivering technology



[linkedin.com/company/bio-strategy](https://www.linkedin.com/company/bio-strategy)

Valid until:
31 Dec 2021
or while stocks last

Ref. Code:
PHCbi_ULT-Freezer_StockSale

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in NZD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.