

LIMITED STOCK - ORDER NOW!

Ultra-Low Temperature TwinGuard Freezers

Takes your valuable sample storage to the next safety level with two independent refrigeration systems. This Dual Cooling System consists of two individual compressors to provide a reliable -86°C Ultra-Low temperature environment. Should an unexpected failure occur in one cooling circuit, the other will maintain at least -70°C.

These freezers represent the ultimate in sample storage and all models feature VIP Plus vacuum insulation, large colour touch screen display, hot gas bypass and 5-fold door gasket to reduce ice build-up, filter-less design – no messy cleaning dirty filters!, automatic vacuum release port, USB port for data download & supplied with 3 adjustable stainless-steel shelves.



MDF-DU302VX-PE

MDF-DU502VX-PE

MDF-DU702VX-PE

Description	Quantity	Code	Price NZ\$	Sale NZ\$
MDF-DU302VX-PE TwinGuard 360L	3	SANYMDF-DU302VX-PE	21,168	15,995
MDF-DU502VX-PE TwinGuard 528L	3	SANYMDF-DU502VX-PE	23,421	18,975
MDF-DU702VX-PE TwinGuard 729L	4	SANYMDF-DU702VX-PE	25,615	19,975

Contact us for more information
and to place your orders

Bio-Strategy Ltd.

T: 0800 34 24 66
E: sales.nz@bio-strategy.com
www.bio-strategy.com
shop.bio-strategy.com



bio-strategy
delivering technology

 [linkedin.com/company/bio-strategy](https://www.linkedin.com/company/bio-strategy)

Valid until:
31 Dec 2021
or while stocks last

Ref. Code:
PHCbi_ULT-Twinguard_StockSale

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in NZD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.