

# Equipment Stock Sale

## Limited Stock Only!



bio-strategy  
delivering technology

### PHCbi Climate Chamber, 294L + Humidity

(Repaired after damage during local transit)

- Control Range:  
+5°C to +50°C (Light OFF); +10°C to +50°C (Light ON)
- Humidity Control Range:  
60 - 90% RH (Light OFF), 55 - 85 % RH (Light ON)
- Forced air circulation
- PID control
- External dimensions/weight 760 x 700 x 835mm (WxDxH), 235 kg
- 6-month warranty



**PHCbi**

Product Code:  
SANYMLR-352H-PE-TD

**1 only**

List Price: 25,459

**Sale Price: 13,500**

### PHCbi CO<sub>2</sub> Chamber, 49L

Personal Type CO<sub>2</sub> Incubator

- Ambient +5°C to 50°C ± 0.25°C
- Thermal Conductivity CO<sub>2</sub> Sensor (0% - 20% ± 0.15%)
- 3 shelves standard
- External dimensions/weight 80 x 548 x 575mm (WxDxH), 49 kg
- 12-month warranty



Product Code:  
SANYMCO-5AC-PE **1 only**

List Price: 7,827

**Sale Price: 5,665**

### PHCbi -30 Chest Freezer, 138L

Outstanding reliability and performance

- Control Range -20°C to -30°C
- Alarms (Hi/Low, Power Failure, Remote Alarm Contacts)
- Includes 2 Baskets, Access Port
- External dimensions/weight 8640 x 687 x 881mm (WxDxH), 52 kg
- 12-month warranty



Product Code:  
SANYMDF-137-PE **2 only**

List Price: 3,913

**Sale Price: 2,886**

### PHCbi Rack 24 Boxes 2ml

Stainless Steel Freezer Rack for 24 x 2 inch boxes  
for MDF-U52V and MDF-U72V

Product Code:  
SANYURO-462

**9 only**

List Price: 250 | **Sale Price: 175**

**PHCbi**



Bio-Strategy Ltd | T: 0800 34 24 663 | E: sales.nz@bio-strategy.com  
NZ\$, excl. GST / local freight, for current stock only, **valid 30 May 22**

# Equipment Stock Sale

## Limited Stock Only!



bio-strategy  
delivering technology

### IKA Rotavisc me-vi complete, Rotational Viscometer package

Accuracy of the ROTAVISC rotational viscometer for both Newtonian and non-Newtonian fluids is +/- 1 % of the measuring range. The reproducibility is +/- 0.2 %.

Includes the following:

- Rotastand
- SP set-2
- PT 100.8
- Guard rail
- USB Cable - USB A to Micro-B, 2m
- 12-month warranty



**IKA®**

Product Code:  
IKAA25000311 **1 only**  
**Sale Price: 4,495**

### IKA RV 10.83 EVAPORATION FLASK (NS 29/32)


The flask is made of high quality borosilicate glass, 500 ml



Product Code:  
IKAA3740400 **1 only**  
List Price: 144  
**Sale Price: 70**

### IKA VX 1 One-hand attachment

- For shaking single, non-fixed vessels of 1 - 250ml
- Dimensions (W x D x H): 130 x 135 x 40mm
- Weight: 160g



Product Code:  
IKAA607200 **1 only**  
**Sale Price: 93**

### IKA S 10 N - 10 G Dispersing Tool

- Suitable for dispersing instrument T 10 basic ULTRA-TURRAX®
- Materials in contact with medium: PTFE, AISI 316L
- Immersion depth: 20 - 100mm
- All Sterilization methods
- Working range: 1 - 100ml
- Stator diameter: 10mm
- Suitable for solvents
- Max. temperature: 180°C



Product Code:  
IKAA3370100 **1 only**  
List Price: 2,052  
**Sale Price: 1,276**

### IKA MS 3.4 Microtiter attachment

For use with a microtiter plate, for the minishakers MS 3 control and MS 3 digital



Product Code:  
IKAA3426400 **1 only**  
List Price: 210  
**Sale Price: 151**

### IKA R 350 Universal Clamp

For clamping flask necks, condenser, etc up to 11 cm diameter



Product Code:  
IKAA1752900 **1 only**  
List Price: 397  
**Sale Price: 138**

PHCbi

rbioPharm

IKA®

Sherwood

grant


TAYLOR-WHARTON

stuart

Bio-Strategy Ltd | T: 0800 34 24 663 | E: sales.nz@bio-strategy.com  
NZ\$, excl. GST / local freight, for current stock only, **valid 30 May 22**

# Equipment Stock Sale

## Limited Stock Only!



### Ridacube Scan Analyser 340/580nm

For Ridacube Test Kits

Run enzymatic and colorimetric assays for the detection of organic acids, sugars or other food components

- RFID technology: Pre-set calibration curve, lot number and expiry, number of tests remaining
- Tablet/PC interface: Android-based app, updates available via internet, USB connection
- Fully automated system but no pipetting device, closed test cartridges, integrated QC self-test
- Weight 2.4 kg, dimensions 16 x 13 x 14.5cm, CE-compliant
- 12-month warranty



Product Code:  
RBIOZRCS0580 **1 only**

List Price: 3,621

**Sale Price: 2,300**

### PCMT Thermoshaker with Cooling for Microtubes and Microplates

- Includes powersupply and 2 x spare rubber belts
- Fast heat-up and cooling times
- Even mixing across the thermoshaker block
- Thermoblock not included
- 6-month warranty



Product Code:  
GRANPCMT **1 only**

List Price: 2,707

**Sale Price: 500**

### Cryogenic Refrigerator LS3000 + 6 Racks

Box/rack system, storing 3000 2ml vials, 100 vial boxes

- Static Holding Time: 106 Days
- Working Time: 66 Days
- Evaporation Rate: 0.76 Liters/Days
- Liquid Nitrogen Capacity: 81 Liters
- Holds: 6 racks with (5)5cm boxes
- Overall Height: 73cm / Overall Diameter: 68cm
- 6-month warranty



Product Code:  
TAYWTAYLS3000 **1 only**

List Price: 10,240

**Sale Price: 2,995**

### Sample Tube for MSB

Normal bore (4mm OD, 3.24 mm ID)

- For use with Sherwood Scientific Magnetic Susceptibility Balance Models MK1 and Auto



Product Code:  
SHEW71086002 **1 only**

List Price: 125

**Sale Price: 80**

### Heating Element for Merit Water Still W4000



Product Code:  
STUAB227-APR **5 only**

**Sale Price: 432**

