

## Axygen® Tips & Tubes

### Safe Guard Your Supply Chain!

**With supply shortages happening worldwide check out our Axygen products to ensure continuance of your testing and research.**

**On Sale**

Axygen® Strip Tubes, Deep Well Plates and Tips for use with Tecan, Qiagen and other Instruments.


0.1ml PCR tubes are specifically designed for use with the Corbett Rotorgene system. There are 4 tubes per strip and 4 caps per strip.

0.2 ml polypropylene thin wall PCR tubes with flat caps are designed to fit standard 96 well blocks and are available in an assortment of colors. Ultra-thin and consistent wall thickness allows precise thermal transfer and even insertion cap ensures a tight fit to prevent evaporation.



RNase-/DNase-free. Nonpyrogenic. Nonsterile. Manufactured with 99.9% pure virgin polypropylene.


Description	Code	UOM	Price	Sale
0.1ml Polypropylene PCR Tube Strips and Caps, Clear, Nonsterile	AXYGPCR-0104-C_U	250 strips = 1,000 tubes	\$148.30	<del>\$120.20</del>
0.1ml Polypropylene PCR Tube Strips and Caps, Clear, Nonsterile	AXYGPCR-0104-C	2,500 strips = 10,000 tubes	\$1,483	<del>\$1,112</del>
0.1mL Low Profile Polypropylene Thin Wall PCR Tube Strips and Real Time Strip Caps, 8 Tubes and Caps/Strip, Clear, Nonsterile	AXYGPCR-0108-LP-RT-C	1,250 strips = 10,000 tubes	\$1,872.20	<del>\$1,404</del>
0.2ml Thin Wall PCR Tubes with Flat Cap, Clear, Nonsterile	AXYGPCR-02-C	1,000 tubes	\$85	<del>\$64</del>



AxyMats™ sealing mats minimize evaporation in a wide range of temperatures in PCR and storage applications. Made of research-grade silicone, they can be autoclaved, and following a bleach wash and ethanol rinse protocol, they can also be reused. For use with 96 square well deep well plates. Compatibility: P-2ML-SQ-C

Description	Code	UOM	Price	Sale
AxyMats™ 96 Square Well Sealing Mat for Deep Well Plates, Nonsterile	AXYGAM-2ML-SQ	50 items	\$384.90	<del>\$281</del>

Axygen® deep well and assay plates are used for sample collection, in vitro growth chambers and long-term storage. Made of polypropylene, they offer excellent chemical resistance and will withstand temperatures down to -80°C. 550 µL deep well plate. ANSI/SLAS 1-2004: Microplates - Footprint Dimensions. Ultra-flat surface ensures proper sealing with heat sealing film.



Description	Code	UOM	Price	Sale
96 Well Clear Round Bottom 550µL Polypropylene Deep Well Plate, 10 per Pack, Nonsterile	AXYGP-96-450R-C	50 items	\$255.80	<del>\$192</del>

# Axygen® Deep Well and Assay Plates

Axygen® deep well and assay plates are used for sample collection, in vitro growth chambers and long-term storage. Made of polypropylene, they offer excellent chemical resistance and will withstand temperatures down to -80°C. 2.2ml deep well plate. ANSI/SLAS 1-2004: Microplates - Footprint Dimensions. Ultra-flat surface ensures proper sealing with heat sealing film.



Description	Code	UOM	Price	Sale
96 Well Clear Round Bottom 2ml Polypropylene Deep Well Plate, 5 per Pack, Nonsterile	AXYGP-2ML-SQ-C	25 items	\$234.90	\$177
96 Well Clear Round Bottom 2ml Polypropylene Deep Well Plate, 5 per Pack, Sterile	AXYGP-2ML-SQ-C-S	25 items	\$348.50	\$262

Axygen® offers a comprehensive line of tips for Tecan®, Aurora Biomed and PerkinElmer® automation systems. These tips are manufactured to stringent specifications under strict process controls and assembled via automation to ensure consistent performance and quality.

Tips are RNase-/DNase-free and nonpyrogenic. Available in nonsterile, sterile, non-filtered and filtered. Clear and Conductive black polypropylene tips are available in 96-well format and is also available in hanging tip rack format designed for Tecan® Genesis Freedom®/Freedom Eva®/Miniprep with LiHa-1, -4, -8 head, Aurora Biomed Versa® with 8-Channel head and PerkinElmer® Janus®/MultiProbe with Variospan-1, 4, 8 head.



Description	Code	UOM	Price	Sale
200µL Sterile Filtered OnCor tip, 96 Tips per rack, 10 Racks per unit, 5 Units/Case.	AXYGQ-OCF-200-R-S	4,800 tips	\$590.20	\$443
1100µL Qiagen BioRobot, Black Conductive 96 Tips per rack, 16 Racks per unit	AXYGQT-1100-CBK-R	1,536 tips	\$368.30	\$277
300uL Qiagen BioRobot, Black Conductive 96 Tips per rack, 24 Racks per unit	AXYGQT-300-CBK-R	2,304 tips	\$614.40	\$461
Axygen Automation Tips for Qiagen® BioRobot® Series, 1000µL, Liquid Level Sensing, Filtered, Sterile 96 Tips/Rack, 16 Racks per unit (max volume 1,100µL)	AXYGQTF-1100-CBK-R-S	1,536 tips	\$737.10	\$553
25µL Pre-sterilised Black PP Conductive Tecan-Style, QIAgility, Liquid Level Sensing Filter Tip, Hanging tip rack, 96 Tips per rack, 24 racks per case (max volume 50µL)	AXYGTTF-50-CBK-HTR-S	2,304 tips	\$744.50	\$559
25µL Pre-Sterilised Black PP Conductive Tecan-Style, QIAgility Liquid Level Sensing Filter Tip, 96 Tips per rack, 2 Racks Per Unit, 36 units per case (max volume 50µL)	AXYGTTF-50-CBK-HTR-S36	6,912 tips	\$1,892.20	\$1,420
175µL Pre-sterilised Black PP Conductive Tecan®, QIAgility Liquid Level Sensing Filter Tip, 96 tips per rack, 24 racks per case (max volume 200µL)	AXYGTTF-200-CBK-HTR-S	2,304 tips	\$440.40	\$331

Contact us for more information and to place your orders

**Bio-Strategy Ltd.**

T: 0800 34 24 66

E: sales.nz@bio-strategy.com

www.bio-strategy.com

shop.bio-strategy.com



**bio-strategy**  
delivering technology

VALID UNTIL:  
31 March 2022

Ref. Code:  
AXYG-Tips



[linkedin.com/company/bio-strategy](https://www.linkedin.com/company/bio-strategy)

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in NZD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.