

Bel-Art No-Wire Test Tube Racks



Submersible! Stackable! Color- Coded! Autoclavable!

A great alternative to wire racks, with price and space saving advantages over other plastic racks.

Dimensions for all racks are 24.6 x 10.4 x 6.4cm high and they assemble in seconds.

- Made of a special blend of autoclavable polypropylene; will sink in a water bath and maintain their stability even in agitated water
- Also perfect for use outside the water bath; meets the demands of daily usage in the lab
- Available in numerous hole sizes, with a choice of three colors for safety and identification
- Features end plates for labeling; numbered and lettered rows for tube identification
- Interlocking feet allow stacking of racks for storage
- Steam autoclavable at 121°C

**Only
\$22!**



Description	Code	Description	Code
13mm Tube White 90 place	BELAF187450000	16mm Tube Blue 60 place	BELAF187470001
16mm Tube White 60 place	BELAF187450001	20mm Tube Blue 40 place	BELAF187470002
20mm Tube White 40 place	BELAF187450002	30mm Tube Blue 21 place	BELAF187470004
13mm Tube Red 90 place	BELAF187460000		

More stock to come. Contact us for pricing.

Description	Code	Description	Code
25mm Tube White 24 place	BELAF187450003	25mm Tube Red 24 place	BELAF187460003
30mm Tube White 21 place	BELAF187450004	30mm Tube Red 21 place	BELAF187460004
16mm Tube Red 60 place	BELAF187460001	13mm Tube Blue 90 place	BELAF187470000
20mm Tube Red 40 place	BELAF187460002	25mm Tube Blue 24 place	BELAF187470003



discount available on current stock only

Contact us for more information
and to place your orders

Bio-Strategy Ltd.

T: 0800 34 24 66
E: sales.nz@bio-strategy.com
www.bio-strategy.com
shop.bio-strategy.com



bio-strategy
delivering technology



[linkedin.com/company/bio-strategy](https://www.linkedin.com/company/bio-strategy)

VALID UNTIL:
31 March 2022

Ref. Code:
BELA-Wireless

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in NZD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.