

Epredia HM325 Manual Rotary Microtome

POWERED BY SHANDON  MICROM  MENZEL-GLÄSER  Richard-Allan Scientific 

Trusted Performance, Added Efficiencies

Ultra-light touch, ergonomic flywheel

- Large wrap-around, detachable waste tray
- Bio-Strategy Microtome service experience



1 unit in stock

- High stability and precision for premium-quality paraffin sections
- Unique section thickness setting with removable knobs allowing operation with both the right and left hands
- Resettable electronic section counter
- Two quick-trim stages (10 µm and 30 µm) independent from the preset fine section thickness
- X/Y fine orientation with reproducible zero positioning
- Manual coarse feed wheel
- Specimen retraction (can be deactivated)
- Long stroke length (64 mm) allows high quality sections of Macro /SuperMega cassettes
- Large, removable section waste tray covers the entire working area

MX35 Ultra™ Low-Profile Blade

All-around excellent low-profile blade for all types of tissue sectioning and increased blade life.

- Soft tissue
- Hard tissue
- Thin sectioning
- Biopsies and recuts

Quantity discounts available for standing orders

Physical Specifications

Dimensions (D x W x H):
19 x 16.5 x 11 in.
(49 x 42 x 28 cm)

Weight: 51 lbs / 23 kg

Horizontal feed range: 28 mm
Section thickness range: 0.5-60 µm
Mechanical trim stages: 10/30 µm
Retraction: 60 µm

Includes

Standard tools, brushed aluminum cover plate, dust cover, and operator guide



Description	Code
HM 325 Manual Rotary Microtome	EPRN902100
With Universal Cassette Clamp and E blade holder	EPRN902100A
With Universal Cassette Clamp and ER blade carrier	EPRN902100ER
MX35 Ultra™ Low-Profile Blade	EPKM3053835

You'll need:



blades,
blade holders,
knives



Cool Cut
and STS



microscope
slides



freezing
sprays

Contact us for more information
and to place your orders

Bio-Strategy Ltd.

T: 0800 34 24 66
E: sales.nz@bio-strategy.com
www.bio-strategy.com
shop.bio-strategy.com



bio-strategy
delivering technology



[linkedin.com/company/bio-strategy](https://www.linkedin.com/company/bio-strategy)

VALID UNTIL:
31 March 2022

Ref. Code:
EPRN-MICRO

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in NZD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.