

Esco Airstream E Class II Cabinets

Free adjustable height stand, UV light and 2 power sockets (valued at \$1,750)

The AC2-E series cabinets comply to the EN12469 standard.

The manual sash automatically activates internal lighting and airflow when set to the marked working height.

ULPA downflow and exhaust filters with 99.999% efficiency, 10X better than HEPA filter (99.99%) improves safety and delivers ISO Class 3 air cleanliness.

- Interlocked UV light with timer – no accidental burns
- Low noise level for a more comfortable working environment
- Also available in 0.6, 0.9, and 1.5m widths
- Ergonomic sloping front combined with tempered glass side panels allows entry of more natural light
- Sentinel Gold microprocessor controller with digital display shows cabinet status
- Multi-piece stainless steel work surface is easy to access and clean
- Night Setback Mode reduces energy consumption by a further 60%
- 3-year warranty

25% Off
+
Free Extras



Description	Code	Price	Sale
1.2m Class II Microbiological Safety Cabinet	ESCOAC2-4E8	\$14,822	\$11,100
1.8m Class II Microbiological Safety Cabinet	ESCOAC2-6E8	\$20,199	\$15,150

Contact us for more information and to place your orders

Bio-Strategy Ltd.

T: 0800 34 24 66

E: sales.nz@bio-strategy.com

www.bio-strategy.com

shop.bio-strategy.com



bio-strategy
delivering technology



[linkedin.com/company/bio-strategy](https://www.linkedin.com/company/bio-strategy)

VALID UNTIL:
31 March 2022

Ref. Code:
ESCO-AC2-E

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in NZD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.