


ESCO Isotherm Ovens

Economical Heating Solutions!

Laboratory ovens are used for high-volume thermal convection applications including annealing, drying, sterilizing, and other lab functions.

**20%
Off!**

- 54, 110 and 240 litre capacities
- Ambient + 7.5°C to 250°C
- Dual display showing set and actual temperature
- Supplied with 2 shelves (additional shelves can be ordered)
- Additional capacities available (32L and 170L)
- Stainless steel interior and shelves
- 4-zone heated air jacket for maximum temperature uniformity
- Forced air circulation with adjustable fan speed
- Air flow slider to adjust exchange rate of air
- Microprocessor PID control with 10-program memory
- Free stacking kit included



Pre-heat Chamber Technology
Ensures stable heating and maximum temperature uniformity in the chamber

Description	Code	List	Price
Isotherm 54L Laboratory Oven	ESCOOFA-54-8	\$3,200	\$2,560
Isotherm 110L Laboratory Oven	ESCOOFA-110-8	\$4,155	\$3,324
Isotherm 240L Laboratory Oven	ESCOOFA-240-8	\$6,878	\$5,495

Contact us for more information
and to place your orders

Bio-Strategy Ltd.

T: 0800 34 24 66
E: sales.nz@bio-strategy.com
www.bio-strategy.com
shop.bio-strategy.com