

Grant Instruments PV-1 Vortex Mixer

In Stock and On Sale!

Ideally sized for personal use

The Grant-bio PV-1 vortex mixer is compact, with a low profile and small footprint, to sit in almost any location.

The PV-1 is so small that it will fit in your hand, yet there is no compromise on performance.

The reliable motor produces regulated and reproducible agitation throughout the speed range.

- Variable speed
- Continuous or touch operation
- Pressure-sensitive cup for tubes up to 20 mm diameter
- For use at 4 to 40°C ambient (80% relative humidity, non-condensing) in a cold room or incubator
- Adjustable speed control: 750 to 3000 rpm
- 'Continuous' or 'touch' operation
- Reliable and extremely quiet motor produces regulated and reproducible agitation throughout the speed range
- Extremely easy to operate – select either 'continuous' or 'touch' operation and turn the knob to adjust the speed
- In 'touch' mode, agitation starts in response to pressure on the pressure sensitive cup.

With variable speed from 750 to 3000 rpm, the unit can be used for gentle through to vigorous mixing and pellet resuspension in 0.5ml, 1.5ml, and 15ml tubes.

Rubber suction pads hold tight to the work surface to prevent the unit from 'walking', and they also absorb vibration which prevents its transmission to the workbench.

SAVE!



Description	Code	Price	Sale
Vortex Mixer PV1 Mini for Tubes to 20mm	GRANPV-1	\$347	\$315

Contact us for more information
and to place your orders

Bio-Strategy Ltd.

T: 0800 34 24 66
E: sales.nz@bio-strategy.com
www.bio-strategy.com
shop.bio-strategy.com