

## Gyrozen Mini Centrifuge

### Personal 12-place Micro Centrifuge

Micro Centrifuge, supplied with autoclavable 12-place fixed angle micro-rotor, snap-fit plastic rotor lid and 12 adaptors for 0.2 ml PCR tubes.

- Max speed: 13,500 rpm; Max force: 12,300 x g
- Microprocessor control with digital display of speed/force & time
- Pulse key and automatic door release
- Quiet in operation with rapid acceleration and deceleration

PCR and cryotube rotor options available.  
Contact us for more details.

**Save!**



Description	Code	Price	Sale
Mini Microcentrifuge, White	GYROGZ-1312	\$1,270	\$999
Mini Microcentrifuge, Blue	GYROGZ-1312BLUE	\$1,329	\$1,100
Mini Microcentrifuge, Pink	GYROGZ-1312PINK		



Contact us for more information  
and to place your orders

**Bio-Strategy Ltd.**

T: 0800 34 24 66

E: sales.nz@bio-strategy.com

www.bio-strategy.com

shop.bio-strategy.com



**bio-strategy**  
delivering technology



[linkedin.com/company/bio-strategy](https://www.linkedin.com/company/bio-strategy)

VALID UNTIL:  
31 March 2022

Ref. Code:  
GRYO-Mini

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in NZD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.