

## 12 MP Microscope Camera Flexacam C3

### New and Improved!

#### Camera simplicity without a PC

Complete your daily imaging tasks efficiently and simply without the need for a PC, including comparison, measurement, annotation, and data sharing.

##### Save time with integrated on-screen display

The Flexacam C3 can be adjusted and operated with the integrated on-screen display. It also offers intuitive measurement, annotation, overlay, and network tools which give more flexibility to enhance documentation.

##### Quickly view, record, and share

- View images directly on your monitor and record videos in full HD with up to 30 frames per second
- Save images to your local network or send as email attachment for easy sharing and storage

##### Enhance your documentation

**Measurement:** Measure multiple sample features in the live image and save the results together with the image

**Overlay:** Facilitate your decision making by comparing the live image directly on a monitor to reference images or customized overlays with a single click

**Annotation:** Highlight specific features and areas of interest on a sample and easily add comments and conclusion to the image

**NEW!**

##### Powerful colour accuracy

Make reliable and accurate decisions for your application tasks with accurate clear images showing the fine details of your sample in true-to-life colour, thanks to the 12 MP CMOS sensor, large dynamic range, and fast signal processing.

##### See colour-precise, crisp images

The Flexacam C3 displays live images fast, with precise colours, and in high resolution for many applications, such as materials science, production and quality control, life science, pathology, clinical, as well as many others.

View images in 4K, or full HD, while capturing fine details of your sample:

- Precise colour reproduction is provided by the sensor's dynamic range and signal processing unit
- Fine details of the sample are revealed thanks to the 12 MP CMOS sensor
- Fast screening of sample with bright and dark areas due to the fast auto exposure



Simplicity, color accuracy,  
and real-time collaboration

Flexacam C3  
12 MP microscope camera


### Collaboration in real time

Collaborate and review samples details efficiently with your team members using the fast display of live images.

### Work efficiently as a team

Several users can view at the same time the live image of the sample which is displayed directly on a monitor. Seeing motion of the sample is the same via the monitor as through the eyepieces. This allows for fast and easy discussion of the sample details being viewed and can save time by making team decisions possible in real time.

- Clear and smooth visualization of sample movements on the screen thanks to the high frame rate of up to 60 fps and high signal transfer speed from sensor to monitor
- Viewing of live and captured images directly on a monitor simultaneously by multiple users when used in stand-alone or workstation mode



### Choose how you connect, save, and share


No need to change how you view, save, and share microscope images. The Flexacam C3 is adaptable to your needs.

### Set it up the way you prefer:

1. HDMI allows the camera to be connected directly to a monitor and used in stand-alone mode without the need for a PC
2. Use the ethernet to work with direct access to your facility's network
3. Save images directly on a USB stick when the camera is used in stand-alone mode
4. To use the camera with a PC, connect them together via the USB port (type-C)
5. Easily connect a mouse, keyboard, or Wi-Fi dongle with the USB ports (type-A)

### Description

- | Description  | Price   |
|--|---------|
| 12 MP Microscope Camera Flexacam C3 - excluding Standalone Kit | \$6,130 |
| 12 MP Microscope Camera Flexacam C3 - including Standalone Kit | \$6,878 |



Contact us for more information  
and to place your orders

### Bio-Strategy Ltd.

T: 0800 34 24 66  
E: sales.nz@bio-strategy.com  
www.bio-strategy.com  
shop.bio-strategy.com



**bio-strategy**  
delivering technology



[linkedin.com/company/bio-strategy](https://www.linkedin.com/company/bio-strategy)

VALID UNTIL:  
31 March 2022

Ref. Code:  
LEIM-FLEX

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in NZD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.