

Miele PG-8583 Washer / Disinfector

Make manual washing a thing of the past!

Package includes:

A102 upper basket with 2 x AK12 inserts,
A150 lower module basket with 2 x A301/5 modules
(3 x 6 = 18 nozzles/module).

Reduce the risk of injury through manual washing and improve safety in your wash area!

- 14 standard programs + 2 freely configurable programs
- Two wash levels configured for washing wide neck and narrow neck glassware.
- Spray arm and spray pressure monitoring ensuring a clean wash
- Underbench or standalone models
- Compact: 600mm x 600mm x 850mm (820mm lid removed)
- Powerful 500L/min pump with integrated heating elements for safe washing cycles up to 93°C
- Passive EcoDry drying function – glassware can be dried
- Stainless steel interior/exterior for durability and cleanliness
- ADP model including pump for AD water and waterproofing kit
- Many other options and accessories to select from
- Dispensing pumps for liquid detergent and liquid neutraliser

Package Deal!



Description	Code	List	Sale
Miele PG-8583 Washer / Disinfector Package 1	MIELPG-8583	\$19,506	\$16,995

Contact us for more information and to place your orders

Bio-Strategy Ltd.

T: 0800 34 24 66

E: sales.nz@bio-strategy.com

www.bio-strategy.com

shop.bio-strategy.com



bio-strategy
delivering technology



[linkedin.com/company/bio-strategy](https://www.linkedin.com/company/bio-strategy)

VALID UNTIL:
31 March 2022

Ref. Code:
MIEL-PG8583

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in NZD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.