

PHCbi Compact -80°C Chest Freezer

In Stock - Great Pricing!

**Ideal Low Temperature Compact Freezer
for vaccine and biological storage.**

**30%
OFF!**

- 84 litres with small footprint: 550 x 685mm
- Adjustable temperature range from -60°C to -80°C
- Filter-less design – reduced maintenance
- Excellent Brand: PHCbi – formerly Panasonic / Sanyo
- Single compressor system, reliable with three year warranty
- Alarms to protect valuable product
- Low noise design; 47dB(A)
- VIP Plus vacuum insulation panel technology
- Improved energy consumption achieving energy savings of approx.. 40%
- Dry “remote” alarm contacts and access port.
- Suitable for operation on 240V 50Hz 10 amp supply



Description	Code	Price	Sale
Compact -80°C Chest Freezer, 84 litres	MDF-C8V1-PE	\$10,288	\$6,975

Contact us for more information
and to place your orders

Bio-Strategy Ltd.

T: 0800 34 24 66

E: sales.nz@bio-strategy.com

www.bio-strategy.com

shop.bio-strategy.com



bio-strategy
delivering technology



[linkedin.com/company/bio-strategy](https://www.linkedin.com/company/bio-strategy)

VALID UNTIL:
31 March 2022

Ref. Code:
SANY-C8V1

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in NZD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.