

Axygen Filter Tips

Quality Axxygen Filter Barrier Tips Designed to Fit most Pipettors

Axygen® Filter Barrier tips are ideal for DNA amplification and microbiological applications where trace contamination can completely invalidate results.

The filter prevents cross-contamination from sample to sample by creating a physical barrier in the tip to inhibit and block aerosols.

Every polyethylene filter, free from DNA inhibitors, is inserted into the tip robotically to ensure freedom from nuclease contamination.

Racked Tips are packaged in sealed-bottom racks with lids to prevent contamination and available in sterilised form.

- RNase-/DNase-free
- Nonpyrogenic
- MaxyMum Recovery Technology for ultra-low retention and high yield
- Sterile with downloadable COA

Choose the MaxyMum Recovery “L-R-S” racked and sterile range for ultra-low retention or “R-S” standard racked and sterile range for greater savings.

For greater reach choose the “MRF-10XT” or “TXLF-10” Extended Length tip.

In Stock!



Description	Pack	Code	Price	Sale
200µL MaxyMum Recovery P200 Rack Sterile [10X96]		AXYGTF-200-L-R-S		
10µL MaxyMum Recovery P2/P10 Rack Sterile [10X96]		AXYGTF-300-L-R-S		
300µL MaxyMum Recovery Universal Fit Rack Sterile [10x96]	960 items	AXYGTF-350-L-R-S	\$116.00	\$86
20µL MaxyMum Recovery P20 Rack Sterile [10X96]		AXYGTF-20-L-R-S		
100µL MaxyMum Recovery Universal Fit Rack Sterile [10X96]		AXYGTF-100-L-R-S		
10µL MaxyMum Recovery Eppendorf Rack Sterile [10X96]		AXYGTF-400-L-R-S		
10µL MaxyMum Recovery Universal Fit Rack [10X96]	960 items	AXYGTXLf-10-L-R-S	\$119.50	\$87
10µL MaxyMum Recovery Extended Length Rack Sterile [10X96]		AXYGMRF-10XT-L-R-S		
1000µL MaxyMum Recovery Universal Fit Rack Sterile [10X100]	1000 items	AXYGTF-1000-L-R-S	\$133.20	\$99
20µL P20 Style Rack Sterile [10X96]	960 items	AXYGTF-20-R-S	\$128.60	\$76
100µL Universal Fit Rack Sterile [10X96]	960 items	AXYGTF-100-R-S	\$117.90	\$74
200µL P200 Style Bulk Non-Sterile	1000 items	AXYGTF-200	\$121.90	\$69

Contact us for more information
and to place your orders

Bio-Strategy Ltd.

T: 0800 34 24 66

E: sales.nz@bio-strategy.com

www.bio-strategy.com

shop.bio-strategy.com



bio-strategy
delivering technology



[linkedin.com/company/bio-strategy](https://www.linkedin.com/company/bio-strategy)

VALID UNTIL:
30 June 2022

Ref. Code:
AXYG-Tips

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in NZD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.