

Corning® Stripettor™ Ultra Pipet Controller and Serological Pipettes

Light and ergonomic and at a lighter price!

In Stock



This light, ergonomic, motorised pipette controller is designed for use with glass and plastic serological pipette from 0.5 to 100ml.

Conveniently positioned switches change operating modes and speeds quickly to handle different liquid volumes and various viscosities.

Up to 8 hours of use when fully charged, with a large LCD display to clearly indicate battery status, pipetting mode, and speed.

Description	Code	Price	Sale
Corning Stripettor Ultra Pipet Controller	CORN4099	\$645.90	\$541.00

Stripette Serological Pipettes

- Color-coding on paper allows easy size identification
- Sterile, individually paper/plastic wrap allows full view
- Unique color-coded magnifying stripe for easier meniscus viewing
- Certified nonpyrogenic and DNase-/RNase-free
- Accuracy within $\pm 2\%$ at full volume
- Lot numbers printed on each pipet wrapper
- Calibrated to deliver (TD-EX 20°C) their contents with blowout



Description	Code	Price	Sale
Stripette Serological Pipette 5ml [200 item]	CORN4487	\$49.70	\$40.70
Stripette Serological Pipette 10ml [200 item]	CORN4488	\$49.70	\$40.70
Stripette Serological Pipette 25ml [200 items]	CORN4489	\$88.60	\$72.70

Contact us for more information
and to place your orders

Bio-Strategy Ltd.

T: 0800 34 24 66

E: sales.nz@bio-strategy.com

www.bio-strategy.com

shop.bio-strategy.com