

Airstream E 1.2m Class II Microbiological Safety Cabinet

In Stock!

**Includes FREE adjustable height stand,
UV light and 2 power sockets (valued at \$1,750)**

The manual sash automatically activates internal lighting and airflow when set to the marked working height.

ULPA downflow and exhaust filters with 99.999% efficiency, 10X better than HEPA filter (99.99%) improves safety and delivers ISO Class 3 air cleanliness.

- Interlocked UV light with timer – no accidental burns
- Low noise level for a more comfortable working environment
- Also available in 0.6, 0.9, 1.5m and 1.8m widths
- Ergonomic sloping front combined with tempered glass side panels allows entry of more natural light
- Sentinel Gold microprocessor controller with digital display shows cabinet status
- Multi-piece stainless steel work surface is easy to access and clean
- Night Setback Mode reduces energy consumption by a further 60%
- 3-year warranty

**20% Off
+
Free Extras**



Description	Code	Price	Sale
1.2m Class II Microbiological Safety Cabinet	ESCOAC2-4E8	\$14,822	\$11,500

Contact us for more information
and to place your orders

Bio-Strategy Ltd.

T: 0800 34 24 66
E: sales.nz@bio-strategy.com
www.bio-strategy.com
shop.bio-strategy.com



 [linkedin.com/company/bio-strategy](https://www.linkedin.com/company/bio-strategy)

VALID UNTIL:
30 June 2022

Ref. Code:
ESCO-AC2E

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in NZD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.