

Esco Compact Isotherm Incubators

Available now !


Esco Isotherm® forced convection laboratory incubators provide a temperature-controlled environment via their forced convection design.

**20%
Off!**

- 54, 110, 170 and 240 litre capacities
- Ambient + 7.5°C to 100°C
- Dual display showing set and actual temperature
- Supplied with 2 shelves (additional shelves can be ordered)
- Sealed glass inner viewing door
- Additional capacities available (32L)
- Stainless steel interior and shelves
- 4-zone heated air jacket for maximum temperature uniformity
- Forced air circulation with adjustable fan speed
- Air flow slider to adjust exchange rate of air
- Microprocessor PID control with 10-program memory



Easy to Clean



Pre-heat Chamber Technology
Ensures stable heating and maximum temperature uniformity in the chamber

Description	Code	List	Price
Isotherm 54L Laboratory Incubator	ESCOIFA-54-8	\$3,306	\$2,645
Isotherm 110L Laboratory Incubator	ESCOIFA-110-8	\$4,332	\$3,465
Isotherm 170L Laboratory Incubator	ESCOIFA-170-8	\$5,693	\$4,555
Isotherm 240L Laboratory Incubator	ESCOIFA-240-8	\$7,037	\$5,630

Contact us for more information
and to place your orders

Bio-Strategy Ltd.

T: 0800 34 24 66

E: sales.nz@bio-strategy.com

www.bio-strategy.com

shop.bio-strategy.com



bio-strategy
delivering technology



[linkedin.com/company/bio-strategy](https://www.linkedin.com/company/bio-strategy)

VALID UNTIL:
30 June 2022

Ref. Code:
ESCO-IFA

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in NZD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.