

Gyrozen Refrigerated Microcentrifuges

Includes FREE 1.5 / 2.0ml Rotor

**Free
Rotor**

Additional rotors also available

- High-speed, quiet and efficient operation
- Automatic rotor ID, imbalance protection & lid lock for enhanced safety
- "Fast Cool" function: 4 °C within 5 minutes!
- Microprocessor control with digital display
- Multiple acceleration & deceleration ramps; pulse key; 100 program memory
- Highly customizable, with a variety of rotors and adaptors available



GYROGZ-1848R-PK
with **FREE** 48-place fixed
angle rotor

- Gyrozen's fastest refrigerated microcentrifuge
- Fixed angle: 18,000 rpm (30,536 xg)
- Swing-out: 13,000 rpm (16,438 xg)
- Can accommodate a range of tube sizes (up to 50ml)
- External Dimension (W x D x H): 380 x 654 x 312mm
- Includes Easy-lock rotor with aluminum lid

Description	Code	Price
1848R Microcentrifuge, Refrigerated with GRE-M-m2.0-48 fixed angle rotor	GYROGZ-1848R-PK	\$8,492



GYROGZ-1730R-PK
with **FREE** 24-place fixed
angle rotor

- Reliable & fast refrigerated microcentrifuge
- Fixed angle: 17,000 (27,237 xg)
- Swing-out: 13,000 (16,343 xg)
- Can accommodate a range of tube sizes (up to 5ml)
- External Dimension (W x D x H): 310 x 692 x 290mm
- Includes Bio-Safe certified Easy-lock rotor for biological containment

Description	Code	Price
1730R Microcentrifuge, Refrigerated with GREB-M-m2.0-24 fixed angle rotor	GYROGZ-1730R-PK	\$6,357

Contact us for more information
and to place your orders

Bio-Strategy Ltd.

T: 0800 34 24 66

E: sales.nz@bio-strategy.com

www.bio-strategy.com

shop.bio-strategy.com



bio-strategy
delivering technology



[linkedin.com/company/bio-strategy](https://www.linkedin.com/company/bio-strategy)

VALID UNTIL:
30 June 2022

Ref. Code:
GYRO-Micro

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in NZD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.