



## Miele Professional PG8504 Washer

### Affordable Washing Solutions

#### Underbench or Standalone Laboratory Washer

With stainless steel interior combining powerful pump with integrated heating elements enabling safe washing cycles up to 70°C.

Miele offers enhanced laboratory safety by removing the manual washing step from your wash process reducing the likelihood of glass breakage that can result in cuts to laboratory staff.

- Door dispensers for powder detergent and neutraliser
- Suitable for operation on 240V 50Hz 13 amp supply
- 6 standard programs + 2 vacant program positions to meet your wash requirements
- Touch on steel controller with 3-line display for easy viewing of washer status
- Compact: 600mm x 600mm x 835mm (W x D x H)
- Stainless steel interior/exterior for durability and cleanliness

Many options and accessories to select from so contact us today for assistance in configuring your washer.

**10% Off!**



Description	Code	Price	Sale
PG8504 Laboratory Washer	MIELG62850422	\$10,620	\$9,550

Contact us for more information  
and to place your orders

#### Bio-Strategy Ltd.

T: 0800 34 24 66  
E: sales.nz@bio-strategy.com  
www.bio-strategy.com  
shop.bio-strategy.com