

## OHAUS Moisture Analysers

### Sleek and Compact Design to Suit Any Workplace

**Intelligent design and effortless operation deliver increased sample throughput, with accurate and reliable results.**

- Flexible Performance, for a wide variety of applications
- Simple to Operate, for faster, smarter and more efficient analysis
- Compact, Easy-to-Clean Design

**20% Off!**

#### MB120: Top Seller

Advanced features to efficiently determine the optimal drying temperature for each sample type, with sufficient storage capacity to allow statistical analysis. Three levels of user management ensure data safety, set up and lock temperature settings and method parameters with administrator access.

#### Key Features:

- SmartGuide™ software: a revolutionary feature that automatically analyses a sample and creates a method, thereby removing the guesswork from the initial setup
- Capacity: 120g; Readability: 0.01%/0.001g
- Temperature Range: 40 to 230°C using halogen lamp
- 4 drying profiles and 7 shut-off criteria
- 110mm TFT touchscreen
- Internal memory capable of storing 100 methods and 1000 measurements



Code	Description	List	Sale
OHAU30246943	MB120 Moisture Analyser	\$5,624	<b>\$4,490</b>

### Accessories



Code	Description	List	Sale
OHAU80252478	Pans, 3 Reusable, 7mm	\$150	<b>\$120</b>
OHAU80252479	Pans, 3 Reusable, 14mm	\$165	<b>\$130</b>
OHAU30585411	Pans, 50, Disposable, Aluminium	\$96	<b>\$76</b>

## OHAUS Moisture Analysers

### MB90:


Geared to increase precision and efficiency in the lab, these Moisture Analysers feature a precisely-controlled halogen heating system which helps to yield fast and repeatable results with readability up to 0.01%. The robust housing and die-cast structure is designed for easy, tool-free cleaning and durability. An intuitive touchscreen icon-driven menu guides users through every step of the operation process.

### Key Features:

- Capacity: 90g; Readability: 0.01%/0.001g
- Temperature Range: 40 to 200°C using halogen lamp
- 2 drying profiles and 3 shut-off criteria
- 110mm TFT touchscreen
- Internal memory capable of storing 2 methods



Code	Description	List	Sale
OHAU30246944	MB90 Moisture Analyser	\$4,721	<b>\$3,775</b>



### MB25:

An ideal combination of results and value, the MB25 is extremely user-friendly. Its compact and easy-to-clean design allows it to fit seamlessly into any workspace.

### Key Features:

- Capacity: 110g; Readability: 0.05%/0.005g
- Temperature Range: 50 to 160°C using halogen lamp
- 1 drying profile and 3 shut-off criteria
- LCD Display

Code	Description	List	Sale
OHAU80252472	MB25 Moisture Analyser	\$3170	<b>\$2,535</b>

Contact us for more information  
and to place your orders

### Bio-Strategy Ltd.

T: 0800 34 24 66

E: [sales.nz@bio-strategy.com](mailto:sales.nz@bio-strategy.com)

[www.bio-strategy.com](http://www.bio-strategy.com)

[shop.bio-strategy.com](http://shop.bio-strategy.com)



**bio-strategy**  
delivering technology



[linkedin.com/company/bio-strategy](https://www.linkedin.com/company/bio-strategy)

VALID UNTIL:  
30 June 2022

Ref. Code:  
OHAU-MOISTURE

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in NZD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.