

Rollex Lec Refrigerators and Freezers

Flexible Storage

12-month warranty as standard

LEC Laboratory Refrigerators and Freezers have been designed with reliability and ease-of-use in mind.

With both under-bench and upright models available, the LEC Laboratory range is SPARK-FREE and suitable for storing combustible supplies safely.

- Lock fitted as standard
- Wire shelves
- Unique anti microbial handle
- Refrigerators operating temperature range +2°C to +10°C
- Freezers operating temperature range -18°C to -20°C

Spark-Free



Description	Code	Price	Sale
LEC Laboratory Refrigerator, Underbench, Solid Door, 158L	ROLLSFSR158	\$3,220	\$3,010
LEC Laboratory Refrigerator, Freestanding, Solid Door, 288L	ROLLSFSR288	\$3,910	\$3,710
LEC Laboratory Freezer, Underbench, Solid Door, 98L	ROLLSFSF98	\$2,800	\$2,660
LEC Laboratory Freezer, Freestanding, Solid Door, 232L	ROLLSFSF232	\$3,850	\$3,650

Contact us for more information
and to place your orders

Bio-Strategy Ltd.

T: 0800 34 24 66
E: sales.nz@bio-strategy.com
www.bio-strategy.com
shop.bio-strategy.com