

Sherwood Scientific Model 926 Mk III Chloride Analyser

Chloride Analyser for the Food and Dairy Industry

Salt measurement (Sodium Chloride) in the Food and Dairy Industry is universal. The ability to read directly, salt concentration and the automatic endpoint detection system has made the 926 the instrument of choice for food manufacturers and analysts throughout the world.

Salt and salinity also represent serious contamination in many industrial processes.

- Robust and versatile coulometric method
- Fast – results in less than 40 seconds! (to 200mg/l)
- Measurement Range: 10 – 999mg/l or 2 – 165mg% Salt
- Ergonomically designed microprocessor-based instrument that is easy to use.
- Optional Active Salt Software (requires PC and balance with RS232) eliminates transcription errors with no need to achieve specific sample and diluent weights.

The sensitivity of the coulometric method enables ppm range measurement in boiler feed water; polymer washes; borehole slurries as well as soil salinity studies. Reading in ppm (mg/l) Chloride, the Select Button offers the immediate conversion of Chloride content into mg% Salt (Sodium Chloride) of the original sample.

Supplied complete with:

- M926 Chloride Analyser
- Printer Cable 9 way RS232
- Buffer (500ml)
- 200mg/l Chloride Standard Solution (100ml)
- 1 Pack of Anodes, 1 Pack of Electrodes
- 2 Marked Beakers
- Silver Electrode Polish
- Stirrer
- Universal Mains Adapter and USB lead A(M) – B(M)
- Operator Manual



Description	Code	Price	Sale
M926 Mk III Chloride Analyser	SHEW92600009	\$11,670	\$9,995

Clinical version (M926S) also available – please contact us for more details.

Contact us for more information
and to place your orders

Bio-Strategy Ltd.

T: 0800 34 24 66

E: sales.nz@bio-strategy.com

www.bio-strategy.com

shop.bio-strategy.com



bio-strategy
delivering technology



[linkedin.com/company/bio-strategy](https://www.linkedin.com/company/bio-strategy)

VALID UNTIL:
30 June 2022

Ref. Code:
SHEW-926Mk-III

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in NZD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.