

Esco Airstream Horizontal Laminar Flow Cabinet, 1.2m

**20%
OFF!**



Including UV Light, PC front cover, 2 electrical outlets and with or without adjustable height stand.

Providing Class 3 clean air work zone by trapping particulates in the ULPA filter.

- Safety ensured: Interlocked UV light with programmable timer
- Glass side panels allowing entry of more “natural” light
- 3-year warranty
- Also available in 0.6, 0.9, 1.2, 1.5, 1.8 and 2.4m widths
- Ergonomic sloping front combined with tempered glass side panels allows entry of more natural light
- ULPA filter with 99.999% efficiency, 10X better than HEPA filter (99.99%) improves safety and delivers ISO Class 3 air cleanliness
- Sentinel Gold microprocessor controller with digital display shows cabinet status
- Single piece spill retaining stainless steel work surface is easy to clean
- Standby Mode reduces energy consumption for better energy utilisation
- Energy efficient low noise ECM motor provides stable self-compensating airflow
- Realtime airflow monitoring
- Isocide™ silver impregnated powder coating inhibits bacterial growth

Description	Code	Price	Sale
With Stand			
Esco Airstream 1.2m Horizontal Laminar Flow Cabinet with Glass Side Panels	ESCOLHG-4AG-F8	\$6,697	
UV Lamp With Safety Interlock and Timer	ESCOUV-30A	\$177	
Polycarbonate Front Cover (Required to Activate UV Lamp)	ESCOFC-LHG/S-4A	\$408	
Power Outlet Inside Cabinet x 2 (\$177 each)	ESCOEO-HE	\$354	\$6,995
Adjustable Height Stand on Lockable Castors	ESCOSTC-4AO	\$1,260	
	Total	\$8,896	
Without Stand			
Esco Airstream 1.2m Horizontal Laminar Flow Cabinet with Glass Side Panels	ESCOLHG-4AG-F8	\$6,697	
UV Lamp With Safety Interlock and Timer	ESCOUV-30A	\$177	
Polycarbonate Front Cover (Required to Activate UV Lamp)	ESCOFC-LHG/S-4A	\$408	\$6,100
Power Outlet Inside Cabinet x 2 (\$177 each)	ESCOEO-HE	\$354	
	Total	\$7,636	

Contact us for more information
and to place your orders

Bio-Strategy Ltd.

T: 0800 34 24 66
E: sales.nz@bio-strategy.com
www.bio-strategy.com
shop.bio-strategy.com

 [linkedin.com/company/bio-strategy](https://www.linkedin.com/company/bio-strategy)

VALID UNTIL:
30 September 2022

Ref. Code:
ESCO-FLOW

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in NZD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.