

PHCbi -30°C Laboratory Upright Freezers

Discount! And FREE TESTO Savaris 2-T3 Wifi data logger with probe.

Available in 3 models to meet your storage requirements.

- 3 capacities: 504L 2-door (upper and lower), 479L 2-door (2 chamber-upper and lower) with independent temperature control of each chamber and 690L single door
- Defrost: manual push button – no temperature spikes within chamber
- 3-year warranty
- Temperature Range: -20 to -30°C
- Filter-less design: no messy filter to clean
- Full alarms + dry remote alarm contacts
- Access ports for probes and sensors
- MDF-MU539-PE and MDF-MU539D-PE supplied with 12 bins and adjustable height shelves
- MDF-U731M-PE supplied with 4 adjustable height shelves

Testo Savaris 2-T3 Digital Wifi Data Logger with T Type Thermocouple included at no charge (normally \$624.00). Testo Savaris 2-T3 can advise end user by email of alarm situation with their freezer. Can log data for up to 3 months. Option to purchase an advanced licence (12/24 or 36 months) for SMS and logging for up to 2 years. Testo Savaris 2-T3 needs an active Wifi network to connect to for email functionality.

20% OFF!



504L

479L

690L

FREE



Description	Code	Price	Sale
PHCbi -30° Upright Freezer with Upper and Lower Doors, 504L incl. Bins	SANYMDF-MU539-PE	\$9,700	\$7,760
PHCbi -30° Upright Freezer 479L Two Door (2 Chamber-Upper and Lower) with Independent Temperature Control of Each Chamber Incl. Bins	SANYMDF-MU539D-PE	\$11,145	\$8,915
PHCbi -30° Upright Freezer 690L Single Door with 4 Adjustable Height Shelves	SANYMDF-U731M-PE	\$12,397	\$9,915
Testo Savaris 2-T3 Wifi Data Logger	TEST572.2033	\$488	FREE
T Type Thermocouple (-200 to +400°C) Incl. Cable	TESTTS419-5MT	\$136	FREE

Contact us for more information and to place your orders

Bio-Strategy Ltd.

T: 0800 34 24 66

E: sales.nz@bio-strategy.com

www.bio-strategy.com

shop.bio-strategy.com