

## Aquastab Algaecide from Scharlau

### Reduce water bath contamination with Aquastab Algaecide!

Aquastab Algaecide is a highly effective blue additive to reduce the contamination from bacteria and fungi in your water bath. Recommended temperature range 5-90°C. Use 5ml/10L of water. 100ml bottle makes up 200L.



**Special  
Price**

Description	Code	Price	Sale
Aquastab Algaecide for Water Bath Use 5ml/10L Water (100ml = 200L)	SCAR0000LZB929	\$112.90	<b>\$95.00</b>

Contact us for more information  
and to place your orders

#### Bio-Strategy Ltd.

T: 0800 34 24 66  
E: sales.nz@bio-strategy.com  
www.bio-strategy.com  
shop.bio-strategy.com

 [linkedin.com/company/bio-strategy](https://www.linkedin.com/company/bio-strategy)

VALID UNTIL:  
30 September 2022

Ref. Code:  
SCAR-AQUA

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in NZD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.