

## DURAN® YOUTILITY®

### The Laboratory Bottle for a Safer Lab!

These Innovative DURAN® YOUTILITY® bottles are safer and easier to handle with their slimmer ergonomic bottle shape and save on space.

- Compare with a Standard Lab bottle to experience the difference!
- Specially shaped gripping zones for easier and safer handling
- New glass thread for faster opening or closing of the YOUTILITY® bottles
- Thread fully compatible with other GL 45 closures and accessories
- Slimmer YOUTILITY® shape optimises the space in autoclaves and refrigerators
- Bottles can be easily stacked horizontally
- Special labelling area compatible with the self-adhesive YOUTILITY® bottle labels
- Easy-to-read ceramic graduation scale indicates nominal volume faster
- Moulded from DURAN® borosilicate 3.3, a pharmacopoeial Type 1 neutral glass for very good chemical and high temperature resistance
- Amber bottle colour meets the USP < 660 > and Ph. Eur. for UV light protection
- Integrated Retrace Codes to allow product traceability via a downloadable certificate of conformance
- Colourful GL 45 Bottle Tags allows for easier identification
- Ergonomic cap shape with comfortable grooves and ridges for easier tightening or removal, especially with smaller sized or gloved hands



**Save  
Today!**

Description	Code	Price	Sale
Bottle 125ml YOUTILITY® Clear W/Cyan Cap and Ring (4 Items)	SCOT218812854	<b>\$45.40</b>	<b>\$37.00</b>
Bottle 250ml YOUTILITY® Clear W/Cyan Cap and Ring (4 Items)	SCOT218813653	<b>\$51.10</b>	<b>\$42.00</b>
Bottle 500ml YOUTILITY® Clear W/Cyan Cap and Ring (4 Items)	SCOT218814452	<b>\$57.10</b>	<b>\$47.00</b>
Bottle 1L YOUTILITY® Clear W/Cyan Cap and Ring (4 Items)	SCOT218815457	<b>\$79.50</b>	<b>\$65.00</b>
Bottle 125ml YOUTILITY® Amber W/Cyan Cap and Ring (4 Items)	SCOT218862859	<b>\$219.90</b>	<b>\$179.00</b>
Bottle 250ml YOUTILITY® Amber W/Cyan Cap and Ring (4 Items)	SCOT218863658	<b>\$262.90</b>	<b>\$214.00</b>
Bottle 500ml YOUTILITY® Amber W/Cyan Cap and Ring (4 Items)	SCOT218864457	<b>\$280.60</b>	<b>\$228.00</b>
Bottle 1L YOUTILITY® Amber W/Cyan Cap and Ring (4 Items)	SCOT218865453	<b>\$395.30</b>	<b>\$322.00</b>

Contact us for more information  
and to place your orders

#### Bio-Strategy Ltd.

T: 0800 34 24 66

E: sales.nz@bio-strategy.com

www.bio-strategy.com

shop.bio-strategy.com

 [linkedin.com/company/bio-strategy](https://www.linkedin.com/company/bio-strategy)

VALID UNTIL:  
30 September 2022

Ref. Code:  
SCOT-YOU

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in NZD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.