

DWK WHEATON® Custom Built Roller Apparatus

Multiple configurations are available so ask for a quote or further information today.

WHEATON® Roller Apparatus, for use with roller bottles, is an established cell culture method and relatively inexpensive way to set up a flexible, scalable biopharmaceutical operation. They are however, heavily in demand, so plan well ahead and order early!

WHEATON®'s Roller Apparatus Team understand that control, monitoring and documentation are vital for every cell culture application. Their state-of-the-art instruments are designed to support industry-standard communication protocols, to enable integration into production facilities. They feature a USB downloadable system event log, for accurate process documentation.

The R2P™ 2.0 controller contains a unique colour touch screen, enabling intuitive operation with gloved hands and offers good distance visibility.



**Order
Today!**

Customise the Apparatus to suit your needs

- New R2P™ 2.0 colour touch screen controller for extremely intuitive operation, or Standard bright LED digital display controller
- Bottom mounted control unit for easier operation by hand and floor-mounted connections, or Top for rigorous cleaning and ceiling-mounted connections
- Production spacing for higher capacity or Modular for easier bottle handling
- Fixed decks with lighter weight for Production use, or Removable decks for greater versatility
- 5 placement areas per deck for use in an upright incubator, or 8 for use in an incubator room for higher capacity
- 1 – 11 decks for optimum capacity
- Optional belt beak alarm and battery backup for the Standard controller units, or optional temperature monitoring, 2 probe alarms and battery backup for the R2P™ 2.0 controller units
- For bottles 108 – 121mm in diameter
- For bottles up to 550mm, however depending on the length of the bottles 1, 2 or even 3 bottles can be used in each placement area
- Standard and CO2 Incubators available to suit most 5 placement models
- Wide range of compatible Wheaton Borosilicate Roller Bottles available

Note: Not all configuration options will work together

Contact us for more information
and to place your orders

Bio-Strategy Ltd.

T: 0800 34 24 66
E: sales.nz@bio-strategy.com
www.bio-strategy.com
shop.bio-strategy.com

Installation, Training IQ/OQ/PQ Validations, Extended Warranty and Service Plans are all available from our Bio-Strategy Service Departments.

 [linkedin.com/company/bio-strategy](https://www.linkedin.com/company/bio-strategy)

Rev. July 2022

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in NZD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.