

Equipment Stock Sale

Limited Stock Only!



bio-strategy
delivering technology

IKA Rotavisc me-vi complete, Rotational Viscometer package

Accuracy of the ROTAVISC rotational viscometer for both Newtonian and non-Newtonian fluids is +/- 1 % of the measuring range. The reproducibility is +/- 0.2 %.

Includes the following:

- Rotastand
- PT 100.8
- USB Cable - USB A to Micro-B, 2m
- SP set-2
- Guard rail
- 12-month warranty



IKA®

Product Code:
IKAA25000311 **1 only**
List Price: \$6,872
Sale Price: \$3,913

IKA VORTEX 3 SHAKER

Vortex shaker suitable for short-time operation (touch function), activated by pressing shaker attachment or continuous operation. Infinitely adjustable speed range: 500 - 2.500 rpm

- Compact design
- Orbital diameter: 4 mm



Product Code:
IKAA3340000 **3 only**
List Price: \$756
Sale Price: \$443

T 25 digital ULTRA-TURRAX®

Stainless steel dispersing elements can be cleaned quickly and easily. Plastic disposable dispersing elements are available in two sizes.

- Digital speed display
- Error code display
- Electronic overload protection
- Quiet operation



Product Code:
IKAA3725000 **1 only**
List Price: \$3,862
Sale Price: \$2,206

IKA S 25 N - 25 G DISPERSING TOOL

Suitable for dispersing instrument T 25 digital
Working range: 50 - 2.000 ml

Stator diameter: 25 mm
Immersion depth: 40 - 165 mm

Suitable for solvents
Max. temperature 180 °C



Product Code:
IKAA1713300 **1 only**
List Price: \$3,101
Sale Price: \$1,743

IKA Color Squid White Magnetic Stirrer

Small magnetic stirrer with glass surface different motives.

- Digital speed display (LED)
- Speed range from 0 - 2.500 rpm
- Max. stirring quantity 1 l
- Recyclable materials



Product Code:
IKAA3671000 **1 only**
List Price: \$670
Sale Price: \$354

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in NZD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.

Equipment Stock Sale

Limited Stock Only!



Ridacube Scan Analyser 340/580nm For Ridacube Test Kits

Run enzymatic and colorimetric assays for the detection of organic acids, sugars or other food components

- RFID technology: Pre-set calibration curve, lot number and expiry, number of tests remaining
- Tablet/PC interface: Android-based app, updates available via internet, USB connection
- Fully automated system but no pipetting device, closed test cartridges, integrated QC self-test
- Weight 2.4 kg, dimensions 16 x 13 x 14.5cm, CE-compliant
- 12-month warranty



Product Code:
RBIOZRCS0580 **1 only**
List Price: \$4,414
Sale Price: \$2,150

PCV-6000 Combined Centrifuge/Vortex Mixer

A variable-speed combined centrifuge/vortex mixer for programmed sequences of spin and mix operations or for independent centrifuging and vortex-mixing of multiple microtubes.

- Compact design
- Centrifugation speed 1000 to 6000rpm
- Optional accessory rotors
- Safety interlock system



Product Code:
GRANPCV-6000 **1 only**
List Price: \$1,493
Sale Price: \$600

SBB Aqua 18 Litre Unstirred Boiling Water Bath

- Adjustable energy regulator provides steady boiling
- Constant level device maintains liquid level
- Robust and reliable design to withstand everyday wear and tear
- Simple, intuitive controls – quick and easy to set temperature
- Non-drip polycarbonate lid included as standard



Product Code:
GRANABBAQUA18PLUS **1 only**
List Price: \$2,591
Sale Price: \$1,100

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in NZD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.



Bio-Strategy Ltd | T: 0800 34 24 663 | E: sales.nz@bio-strategy.com
valid until 23 December 22 or while stocks last