

Polypropylene Tips RNase-/DNase-Free

Universal Tips are designed to fit a wide variety of single and multi-channel pipettors, compatible with most popular brands of pipettors.


These polypropylene tips are nonpyrogenic and RNase-/DNase-free.

Our Axygen tips are available in Extended Length, Maxymum Recovery (low retention), Eppendorf-style, Wide-bore, Sterile and Non-sterile, Filtered and Non-filtered, Racked, Refillable Rack System and Bulk to suit your needs.

Flexible walls and a series of internal sealing rings ensure a secure fit with less force required to load and eject.

1000 μ L Beveled, Blue, Nonsterile

| Description | Code | Sale |
|---|--------------|----------------|
| 1000 μ L Blue Universal Fit [1,000 ITEMS] | AXYGT-1000-B | \$21.00 |



200 μ L, Beveled, Yellow, Nonsterile

| Description | Code | Sale |
|--|-------------|----------------|
| 200 μ L Yellow Universal Fit [1,000 ITEMS] | AXYGT-200-Y | \$18.50 |

VALID UNTIL:
31 Mar 2023
Ref. Code:
AXYG-Tip



Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in NZD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.