





Energy-Efficient
DC ECM Motor



Sentinel™
Gold Microprocessor Controller



Ergonomic
Angled Front

Airstream E Class II Cabinets Type A2

Includes FREE adjustable height stand, UV light and 2 power sockets

Description	Code	Price	Sale
1.2m Cabinet	ESCOAC2-4E8	\$17,730	\$12,995
1.8m Cabinet	ESCOAC2-6E8	\$24,195	\$17,750

Providing operator, product and environment protection, The manual sash automatically activates internal lighting and airflow when set to the marked working height.

- Interlocked UV light with timer – no accidental burns
- Two internal power sockets
- Low noise level (only 53.5dB(A) and 55.7dB(A) respectively)
- Also available in 0.6, 0.9, and 1.5m widths
- Ergonomic sloping front combined with tempered glass side panels allows entry of more natural light
- ULPA downflow and exhaust filters with 99.999% efficiency, 10X better than HEPA filter (99.99%) improves safety and delivers ISO Class 3 air cleanliness
- Sentinel Gold microprocessor controller with digital display shows cabinet status
- Multi-piece stainless steel work surface is easy to access and clean
- Night Setback Mode reduces energy consumption by a further 60%
- Adjustable height stand has lockable castors
- 3-year warranty



VALID UNTIL:
28 Feb 2023
Ref. Code:
ESCO-ClassII

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in NZD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.

