

406 Low Speed Centrifuge

Supplied complete with a fixed-angle rotor including sleeves for 15 ml tubes

The compact and robust 406 centrifuge is ideal for routine work with small sample sizes.

Description	Code	Price	Sale
406 Centrifuge With 6 X 15ml Fixed-Angle Rotor	GYROGZ-0406	\$2,350	\$1,995

- Max. speed: 4,000 rpm / 2,075 x g with 6 x 15ml
- Smart design to minimise heat and noise transmission
- Program memory for up to 10 programs
- Timed runs of up to 99 minutes or continuous operation
- Easy-to-read LCD display



Fixed Angle Rotor

VALID UNTIL:
31 Jan 2023
Ref. Code:
GYRO-406



Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in NZD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.