




The Sterile Sampling Bag That Meets All Your Needs

A secure, contaminant free pliant container that ensures dependable analysis results. You can count on TWIRLEM® products as a proven, economical and efficient way to collect, contain and carry samples.

Description	Code	Sale
140X229mm 76um 650ml - Clear [500 items]	LABPEPL-5590	\$126
178X305mm76um 1650ml - Clear [250 items]	LABPEPL-7012	\$90

Their built-in secure closure tabs are ideal for the shipping and storage of all forms of sample material.



Step1
Tear off the top of the bag along perforations



Step 2
Use pull-tabs to pull open the bag



Step3
Place your sample in to the bag



Step 4
Hold bag by wire ends and twirl to close



Step 5
Fold in ends to provide air-tight and leak-proof closure



products shown with white write on strip are EPR version which is also available.

VALID UNTIL:
31 Mar 2023

Ref. Code:
LABP-EPL

