



EcoDry via Patented
AutoOpen Function



Optimum Wash
Pressure



Crevice-free
Chamber



Spray
Arm Sensing

PG8583 Washer/Disinfector Package

Package includes:

A102 upper basket with 2 x AK12/1 inserts, A150 lower module basket with 2 x A301/5 modules (3 x 6 = 18 nozzles/module).

| Description | Code | Price | Sale |
|---|--------------|-----------------|-----------------|
| PG8583 ADP Underbench / Stand Alone Laboratory Washer / Disinfector | MIEL62858306 | \$20,916 | \$18,500 |

- 14 standard programs + 2 freely configurable programs
- Two wash levels configured for washing wide neck and narrow neck glassware
- EcoDry passive drying
- **Free** detergent and neutraliser and **free** commissioning (once customers local plumber and electrician has connected unit to required services)
- Spray arm and spray pressure monitoring ensuring a clean wash
- Compact: 600mm x 600mm x 850mm (820mm with lid removed)
- Underbench or standalone model
- Powerful 500L/min pump with integrated heating elements for safe washing cycles up to 93°C
- Passive EcoDry drying function – glassware can be dried
- Stainless steel interior/exterior for durability and cleanliness
- ADP model including pump for AD water and waterproofing kit
- Many other options and accessories to select from
- Dispensing pumps for liquid detergent and liquid neutraliser



VALID UNTIL:
28 Feb 2023

Ref. Code:
MIEL-PG8583



Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in NZD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.